

PURELAB Quest



Operator Manual MANU41196 VERSION 02 2/20 Product No. PQDIUVM1 Product No. PQDIXXM1

WATER PURIFICATION

QUICK START GUIDE

WELCOME

Thank you for purchasing the PURELAB® QUEST. Please read the following Health and Safety notice before proceeding.



ENGLISH



SYSTEM CONFIGURATION



ELGA 🕡 VEOLIA

SET - UP / COMMISSIONING

POSITIONING PURELAB QUEST

The PURELAB QUEST unit should be installed on a flat, level surface, in a clean dry environment 15 - 25 °C. The unit can also be wall mounted using the LA735 wall mounting kit against a vertical wall capable of supporting the combined weight of the unit and mounting kit.



b) Place the product on a flat level surface, remove from the packaging. (Fig.2) a) Use lifting aids to move the product into position. (Fig.1) Health and Safety Executive, INDG143 (REV 3), Open Government License,



RECOMMENDED DRAIN AND OVERFLOW CONNECTIONS

The drain and overflow tubes will require a drain or a sink capable of handling at least 1.5L / Min with a gravity fall below the level of the unit. Any connections directed to drain or sink should have an air-break device fitted.

a) Use 6 X flow bends from the basic installation kit LA762.





SET - UP / COMMISSIONING

4 CONNECT THE WATER

RO/DI/UV - Potable Water Supply RO/DI - Potable Water Supply See Operator's Manual - Feedwater Specifications for more information.



- a) Remove the transit plug from the water inlet connection. (Fig.9)
- b) Using the tube supplied, firmly push one end securely into the water inlet connection. (Fig.10)



(Fig.11)

a) Remove the transit plug from the Overflow connection. (Fig.11)



b) Using the tube supplied, firmly push one end securely into the Overflow connection. (Fig.12)



a) Remove the transit plug from the Drain connection. (Fig.13)



 b) Using the tube supplied, firmly push one end securely into the Drain connection. (Fig.14)

Connect the other end of the drain and overflow tubes to a sink or suitable drain capable of handling at least 90L/hr. The drain point should have a gravity fall below the level of the unit and any connections directed to drain should have an air break device fitted. (See Quick Start Guide - Step 3 for details)

C. Turn on the Feedwater supply, open 15mm emergency shut-off valve and regulator the pressure to 4 bars.

ADDITIONAL PURIFIED WATER OUTLET CONNECTIONS

Ports 4,5,6 and 7 do not need to be connected for a standard installation (See Operator's manual - Section 6 for details).



SET - UP / COMMISSIONING



g) Please wait for PURELAB QUEST 7L internal reservoir to fill and reach full purity before dispensing any purified water. If an additional external reservoir is also connected this will start to fill once the internal reservoir is full.



CON	TENT	
	Quick Start	
1.	Health and Safety Notes	Page 1
1.1	Environnent	Page 1
1.2	Electricity	Page 1
1.3	Ultraviolet Light	Page 1
1.4	Residual Biostat	Page 1
1.5	Personal Protective Equipment (PPE)	Page 1
1.6	Control of Substance Hazardous to Health (COSHH)	Page 1
2.	Introduction	Page 2
2.1	Health and Safety	Page 2
2.2	Product Models	Page 2
2.3	Use of this Manual	Page 2
2.4	Installation	Page 2
2.5	Environment	Page 2
2.6	Commissioning	Page 2
2.7	Customer Support	Page 2
3.	Guide to Your PURELAB [®] OUEST	Page 3
4.	Product Consumables and Accessories	Page 4
5.	Register Your Product	Page 5
6.	Additional Purified Water Outlet Connections	Page 6 - 8
6.1	Port 4: Type 2 - Water Outlet Connection for an application	Page 6
6.2	Port 5: Type 3 - Water Outlet Connection for an external reservoir	Page 7
6.3	Port 6 and 7 - Water Outlet Connection for a Free-standing Dispenser	Page 8
7.	Quick Reference Guide	Page 9 - 10
7.1	Manual Dispense	Page 9
7.2	Auto Volume Dispense	Page 9
7.3	Deactivating - Auto Volume Dispense	Page 10
7.4	ON - OFF, Menu and Accept	Page 10
8.	Maintenance	Page 11 - 21
8.1	Replacing Ultraviolet Lamp LC210	Page 11 - 13
8.2	Replacing DI Purification Pack LC292	Page 14
8.3	Cleaning the Inlet Filter	Page 15 - 16
8.4	Replacing the Reverse Osmosis Module LC217	Page 17
8.5	Replacing the Vent Filter LC291 and Point-of-Use Filter LC134 / LC145 / LC197	Page 17 - 18
8.6	Sanitization Procedure - LC209 Sanitization Pack	Page 19 - 20
8.7	Sanitization Procedure - LC209-M2 Sanitization Pack	Page 20 - 21
9.	Troubleshooting	Page 22 - 23
10.	Customize Operation	Page 24
10.1	Set Alarm Points	Page 24
10.2	Advanced Data Logging	Page 24
10.3	Software Update	Page 24
11.	Technical Specifications	Page 25 - 26
12.	Warranty / Conditions of sale	Page 27
Copyright N	ote	
The	information contained in this document is the property of VWS (UK) Ltd, trading as ELGA Veolia®, and is supplied without li	ability for errors or omissions.
	No part of this document may be reproduced or used except as authorized by contract or other written permission fr	rom VWS (UK) Ltd.
	The copyright and an restrictions on reproducing and use apply to all media in which this information may b VWS (UK) Ltd. Pursue a policy of continual product improvement and reserve the right to alter without notice the sp	ecification, design,

price or conditions of supply of any product or service.

@ VWS (UK) Ltd. 2020 - All rights reserved. ELGA $^{\! \otimes}$ and PURELAB $^{\! \otimes}$ are registered trademarks of VWS (UK) Ltd.

Copyright Note

Effersan^{TM} is a registered trademark of Activon Inc.

I. HEALTH AND SAFETY NOTICES

Tape Measure

SanitizationTablets



Adjustable Spanner



INTRODUCTION

2.1 Health & Safety

Please ensure you read the Health & Safety notes in the front of this Operator Manual.

2.2 Product Range

This Operator manual has been prepared for PURELAB® Quest product models:

- PURELAB® Quest RO/DI/UV Product No. PQDIUVM1
- PURELAB® Quest RO/DI Product No. PQDIXXM1

PURELAB[®] Quest converts tap water directly to Type 1, 2 and 3 lab water. The unit is specifically designed for low volume applications and the intermittent demand of the modern laboratory.

2.3 Use of this Manual

This manual guides you through the operation and maintenance of PURELAB® QUEST allowing you to obtain a guaranteed supply of purified water to meet your requirements.

2.4 Installation

The Quick Start Guide shown in the front of this manual shows you how to install and commission the <u>PURELAB®</u> QUEST.

2.5 Environment

The PURELAB[®] Quest unit should be installed on a flat, level surface, in a clean, dry environment. The unit can also be wall mounted using the purpose designed wall mounting kit against a vertical wall capable of supporting the combined weight of the unit and mounting kit LA735. See INST39088 in-box wall mounting instructions for details.



WARNING!

NOT TO BE USED IN FUME EXTRACTION ENCLOSURES.

2.6 Commissioning

<u>PURELAB®</u> Quest is supplied in a preset commissioning mode which must be completed before purified water can be dispensed, see quick start to complete this procedure.

2.7 Customer Support

If you need help with your product, please call your local ELGA Veolia[®] representative. For the address of the nearest ELGA Veolia[®] Sales and Service office visit our website.

www.elgalabwater.com or contact ELGA Veolia at: Email: techsupport@elgalabwater.com



3. GUIDE TO YOUR PURELAB QUEST



4. PRODUCT CONSUMABLES AND ACCESSORIES

Consumables				
Part No.	Description	Typical Service Life*	Max. Shelf Life	
LC292	Purification Pack (DI)	6 Months	2 Years	
LC209	Sanitization Pack	12 Months	2 Years	
LC209-M2	Tablet Sanitization Pack	12 Months	2 Years	
LA728	Bypass Pack (Removed during commisioning)	N/A	N/A	
LC210	UltraViolet Lamp (UV)	12 Months	2 Years	
LC291	Vent Filter (VF)	1 Year	2 Years	
LC217	Reverse Osmosis Module (RO Service Module)	Typical Life 2 Years	2 Years	
LC134	Point-Of-Use 0.2µm Mirco Filter (POU)	3 Months	2 Years	
LC145	Point-Of-Use 0.2µm Mirco Filter (POU)	3 Months	2 Years	
LC197	Point-Of-Use Bio Filter (POU)	3 Months	2 Years	
SKU:40020	Effersan [®] Sanitization Tablet	1 every 6 months	8 Months	
Accessories				
LA821	External Pre-treatment Filter Kit	6 Months	2 Years	
LA735	Wall Mounting Kit	N/A	N/A	
LA512	Feedwater Pressure Regulator	N/A	N/A	
LA757	External 15L Reservoir	N/A	N/A	
LA758	External 30L Reservoir	N/A	N/A	
LA648	External Boost Pump (Feedwater less than <2 Bar)	N/A	N/A	

* Consumable Life is an estimate only, and will depend on the application and feedwater quality.





5. REGISTER YOUR PRODUCT

5. Registering a Product

Taking the time now to register will mean that we can provide a better service to you in the future. We can contact you about product information and service updates.

Why Register your product?

- <u>Validated your product warranty</u>
- Proof of product registration
- <u>Receive software and service updates</u>

How can I quickly register?

The model number and serial number can be found on the ratings plate, see below:



* Model number

** Serial number

Enter the model and serial number into the online form at: https://www.elgalabwater.com/support/register-a-product

6. ADDITIONAL PURIFIED WATER OUTLET CONNECTIONS

The purified water outlet connections are located at the rear of PURELAB[®] QUEST, these are used for purified connections to other applications or products.

Caution! Port 4: Type 2 Outlet should not be connected to a pressurized application, this may cause back-flow into PURELAB[®]Quest.

6.1 <u>PURELAB Quest Port 4: Type 2 - Water Outlet Connection for an application</u>

STEP 1 - Removing the transit plug and connecting the tube

1. Remove the transit plug from water outlet connection Port 4: Type 2 Outlet (Fig.1).

2. Insert an 8mm plastic tube into Port 4: Type 2 Outlet (Fig.2).



STEP 2 - Opening the Purified Water Outlet Valve

- 1. Undo the quarter turn fastner to the access panel (Fig.3).
- 2. Connect the 8mm tube to an application or product that requires purified water.
- 3. Open the manual isolation valve (Turn 90° Anti-clockwise) (Fig.4).



Open

Ţ



6. ADDITIONAL PURIFIED WATER OUTLET CONNECTIONS

6.2 PURELAB Quest Port 5: Type 3 - Water Outlet Connection for an external reservoir

STEP 1 - Removing the transit plug and connecting the tube

- 1. Remove the transit plug from water outlet connection Port 5: Type 3 Outlet (Fig.1).
- 2. Insert an 8mm plastic tube into Port 5: Type 3 Outlet (Fig.2).



STEP 2 - Connecting water and power to the reservoir

- 1. Remove the transit plug from reservoir Port 1: Inlet (Fig.3).
- 2. Insert the 8mm plastic tube from PURELAB Quest Port 5: Outlet into reservoir Port 1: Inlet (Fig.3).
- 3. Insert an 8mm plastic tube into the overflow and connect to drain (Fig.4).



4. Connect the supplied comms lead from PURELAB Quest into J2 reservoir comms connection (Fig.5).





6. ADDITIONAL PURIFIED WATER OUTLET CONNECTIONS

6.3 Port 6 and 7 - Water connections for Free-Standing Dispenser

STEP 1 - Disconnect the power

1. Disconnect the power supply before removing the metal loop (Fig.1).



WARNING! MAKE SURE THAT THE POWER SUPPLY CANNOT COME INTO CONTACT WITH WATER.

STEP 2 - Removing the metal loop

1. Wait a couple of minutes before removing the metal loop.



- 2. Remove the metal loop from Port: 6 and Port: 7 using an 8mm tubing key (Fig.2).
- 3. Follow the Free-standing Dispenser Quick Start Guide to complete this installation.



7. QUICK REFERENCE GUIDE



Press dispense at anytime, release to stop.

7.2 Auto Volume Dispense 🖁



Press Auto-Volume at anytime.



Press to scroll up and down, press dispense to confirm.



7. QUICK REFERENCE GUIDE

7.3 Deactivating - Auto Volume Dispense 🗴



Press Auto-Volume dispense at anytime to exit.

7.4 ON - OFF, Menu and Accept



Deactivate Auto-Volume dispense, press once.



Press and Hold for 2 seconds to enter the Main Menu. Scoll to select a feature in the Main Menu.



Maintenance

Note:

Any maintenance work not included in this operator manual should be carried out by an approved Service engineer or technician, please contact: https://www.elgalabwater.com/support/product-support



Disposing of used consumable items should be in accordance with local statutory regulations.

8.1 Replacing Ultraviolet Lamp LC210 - PQDIUVM1

The Ultraviolet (UV) Lamp should be replaced when;

- UV Lamp failure occurs.
- UV Lamp efficiency decreases affecting water quality.

If either situation occurs, a message will be displayed on screen asking you to replace the UV Lamp.



IF BROKEN MAY CAUSE REDNESS OR IRRITATION AS A RESULT UPON CONTACT WITH SKIN OR EYES, PLEASE SEEK MEDICAL ATTENTION.



STEP 1 - Switching OFF the power supply

1. ISOLATING the power supply, follow the instructions below (Fig.1).





STEP 2 - Removing LC210 UV Lamp

- 1. Open the left side door (Fig.2).
- 2. Disconnect the UV lamp wiring connector (Fig.3).





(Fig.2)

(Fig.3)

STEP 3 - Removing LC210 UV Lamp continued



- 1. Carefully remove the UV lamp from it's chamber (Fig.4).
- 2. Discard used UV lamp in accordance with local regulations (Fig.5).





STEP 4 - Replace LC210 UV Lamp

1. Unpack the new LC210 lamp from it's product packaging (Fig.6).





- 2. Wipe the glass with the alcohol wipe
- 3. Slide the new LC210 UV lamp into the UV chamber (Fig.7).
- 4. Connect the UV lamp connector, making sure it is aligned correctly with the pins (Fig.8).
 - (Use the black arrows as a guide) Alcohol Wipe New LC210 Lamp PURELAB



STEP 5 - Switching ON the Power Supply

1. Connect the power supply. (Fig.9).





STEP 6 - Resetting the LC210 UV Lamp reminder

- 1. Hold and press accept \leftarrow to enter the main menu. (Fig.10).
- 2. Scroll \wedge + \vee and select 'Reset Reminders', press accept \leftarrow to confirm (Fig.11).
- 3. Scroll and select the UV filter, press accept \leftarrow to confirm (Fig.12).
- 4. The reset reminder is now set, normal operation will resume (Fig.13).



(Fig.10)





(Fig.12)



(Fig.13)

PURELAB® Quest Operator Manual MANU41196 Version 02 2/20 Page 13



(Fig.4)

8. MAINTENANCE

8.2 Replacing LC292 DI Purification Pack

Step 1 - Removing the used LC292 Purification Pack

- 1. Press and hold accept ← for 2 seconds to enter the Main Menu (Fig.1).
- 2. Scroll $\wedge + \vee$ and press accept \leftarrow 'Changing Purification Pack'
- 3. Open the right hand door (Fig.2).
- 4. When instructed release the LC292 Purification Pack by 'Pull back' and 'Push down' on the handle (Fig.3).
- 5. Dispose of in accordance with local regulations (Fig.4).



(Fig.3)

(Fig.1)

(Fig.2)

Step 2 - Unpacking the New LC292 Purification Pack

1. Unpack the new Purification Pack and remove the transit plugs (Fig.5).



Step 3 - Installing the New LC292 Purification Pack

- 1. When instructed, insert the base of the new Purification Pack first, then push down and inwards (Fig.6).
- 2. Follow the on screen instructions to complete the procedure.



(Fig.6)



8.3 Cleaning the Inlet Feedwater Filter

PURELAB® QUEST is fitted with an inlet filter to protect it from water particles, which could affect the peformance of the product. It is recommended that the inlet filter assembly is periodically cleaned to remove any particles or debris that can cause a blockaged.

STEP 1 - Isolate the power supply

1. Isolate the power supply to the product (Fig.1).



STEP 2 - Turn OFF Feedwater Supply

STEP 3 - Dismantle the Inlet Filter Assembly

- 1. Remove the feedwater plastic tube, using an 8mm tubing key (Fig.2).
- 2. Unscrew the inlet filter nut with an adjustable spanner (Fig.3).





STEP 4 - Cleaning the Inlet Feedwater Filter

- 1. Check the mesh filter for damage and clean with purfied water (Fig.4).
- 2. Insert the black rubber washer and mesh filter into its original position
- (Ensure that the mesh filter dome is facing outwards)
- 3. Re-Assemble the inlet feedwater filter assembly (Fig.5).



STEP 5 - Reconnecting the Feedwater Supply and Power Supply

- 1. Reconnect the tube and open the valve, restoring the water connection.
- 2. Reconnect the power supply to the product (Fig.6).
- 3. Switch ON to resume normal operation.





8.4 Replacing LC217 Reverse Osmosis Module

PURELAB® QUEST is fitted with a Reverse Osmosis (RO) module which is replaced during a product service or if permeate water purity or flow rate is not adequate and does not meet predicted performance levels. Please contact product support for more information: https://www.elgalabwater.com/support/product-support

8.5 Replacing Vent Filter LC291 and Point-of-Use Filter LC134 / LC145 / LC197

Step 1 - Removing the used LC291 Vent Filter

- 1. Unscrew the used LC291 vent filter from the reservoir (Fig.1).
- 2. Dispose of the used LC291 vent filter in accordance with local regulations (Fig.2).







- 1. Unpack the new LC291 Vent filter from its outer packaging, write the installation date (Fig.3).
- 2. Screw the new LC291 Vent filter into the reservoir (Fig.4).



(Fig.3)

(Fig.4)

Step 3 - Reset the Vent Filter Reminder

- 1. Hold and press accept \leftarrow to enter the main menu. (Fig.5).
- 2. Scroll \wedge + \vee and select 'Reset Reminders', press accept \leftarrow to confirm (Fig.6).
- 3. Scroll and select the VF, press accept ← to confirm (Fig.7).
- 4. The reset reminder is now complete, normal operation will resume (Fig.8).





Replacing the Point-Of-Use Filter LC134 / LC145 / LC197

Step 1 - Removing the Dispense Tip

- 1. Unscrew the dispense tip from underneath the display (Fig.1).
- 2. Do not dispose of the dispense tip, this part is required for the sanitization procedure (Fig.2).



(Fig.1)

(Fig.2)

Step 2 - Unpacking the LC134 / LC145 / LC197 Point-Of-Use Filter and Installing

- 1. Unpack the new POU filter from it's outer and inner packaging, remove the cap (Fig.3).
- 2. Coupler is not required for PURELAB® QUEST, screw in the Point-of-Use filter (Fig.4).



(Fig.4)

- 3. Slightly open both bleed valves, before dispensing any water to release any trapped air (Fig.5).
- 4. Press the dispense button, drawing water through and close both bleed valves (Fig.6).
- 5. Installation is completed, normal operation can resume.





The unit is sanitized to maintain or recover microbiological purity of the product water, sanitization is recommended every 12 months.

- It should be carried out if the product has been switched OFF for a prolonged period.
- If the unit has been operated under adverse conditions such as in high temperatures.
- If the bacteria counts are outside the recommend levels See section 11 Technical Specifications.



BEFORE STARTING THIS PROCEDURE LABEL PURELAB® QUEST WITH THE ATTACHED WARNING SIGN.

8.6 Sanitization Procedure - LC209 Sanitization Pack

Step 1 - Unpacking the LC209 Sanitization Pack

- 1. Remove the Sanitization Pack LC209 from the outer packaging.
- 2. Remove the tranisit plugs from either end of the pack and dispose (Fig.1).



Step 2 - Remove the used Purification Pack LC292 and Install the LC209 Pack

- 1. Press and hold accept \leftarrow for 2 seconds to enter the Main Menu (Fig.2).
- 2. Scroll \wedge + \vee and select 'Sanitization' \leftarrow in the Main Menu (Fig.2).
- 3. Open the right hand door to the Purification Pack housing (Fig.3).
- 4. Release the pack 'Pull back' and 'Push down' on the handle (Fig.4).
- 5. Dispose of the used Purification Pack in accordance with local regulations (Fig.5).



PURELAB® Quest Operator Manual MANU41196 Version 02 2/20 Page 19



Step 3 - Installing LC209 Sanitization Pack

- 1. Insert the base of the Sanitization Pack LC209 first, then clip into position (Fig.6).
- 2. Follow the on screen instructions to complete this procedure.









Locked

8.7 Sanitization Procedure - LC209-M2 Sanitization Pack (USA ONLY)

Step 1 - Unpacking the LC209-M2 Sanitization Pack

- 1. Remove the Sanitization Pack LC209-M2 from the outer packaging.
- 2. Remove the tranisit plugs from either end of the pack and dispose (Fig.1).



Step 2 - Removing the base from LC209-M2

1. Unscrew the pack base, using the red molded hand tool (Fig.2).





Step 3 - Insert a Sanitization Tablet into LC209-M2

- 1. Insert 1 X Effersan Sanitization tablet into the chamber (Fig.3).
- 2. Re-assemble the LC209-M2 pack (Fig.4).
- 3. Hand tighten using the provided moulded tool (Fig.5).
- 4. Ensure there are no gaps in the pack base (Fig.6).



Step 4 - Remove the used Purification Pack LC292 and Install the LC209-M2 Pack

- 1. Press and hold accept \leftarrow for 2 seconds to enter the Main Menu (Fig.7).
- 2. Scroll $\wedge + \vee$ and select 'Sanitization' \leftarrow in the Main Menu (Fig.7).
- 3. Open the right hand door to the Purification Pack housing (Fig.8).
- 4. Release the pack 'Pull back' and 'Push down' on the handle (Fig.9).
- 5. Dispose of the used Purification Pack in accordance with local regulations (Fig.10).



Step 5 - Installing LC209-M2 Sanitization Pack

1. Insert the base of the Sanitization Pack LC209-M2 first, then clip into position (Fig.11).

2. Follow the on screen instructions to complete this procedure.



(Fig.11)



9. TROUBLESHOOTING

This section covers some of the issues that could occur with your PURELAB[®] QUEST, and provides the troubleshooting steps to help you resolve the problem.

For product support, please contact: https://www.elgalabwater.com/support/product-support



WARNING!

ALWAYS ENSURE THAT THE POWER SUPPLY IS ISOLATED, DURING ANY MAINTENANCE WORK.

Issues	Action	
	1. Press any button to wake from sleep mode.	
	or	
Blank Display Screen	2. Check the power supply is connected to the product, ensuring that all connections are securly fitted. (Turn the power OFF/ON)	
	or	
	3. Check the plug socket fuse, replace if nessary.	
	1. Check the product is turned ON.	
Unable to Enter the Features in the	or	
Main Menu	 Make sure Auto-Volume is disabled. (See section 7.3 - Quick Reference Guide) 	
Message Displayed - Purification Pack Change Reminder	1. Replace the Purification Pack. (See section 8.2 Replacing Purification Pack LC292)	
Incorrect Volume of Water	1. Select Auto-Volume Calibration in the main menu.	
	1. Do not exceed 10L of dispened water a day. * Available Volume of Type 1 and 2 water, combined; Increased use will reduce purification pack life.	
Durification Deals Life	or	
Purincation Pack Life	2. Test the feedwater supply, using a water testing kit.	
	Contact product support for advice: https://www.elgalabwater.com/support/product-support	
	(See section 11 Technical Specifications - Feedwater - Conductivity)	
	or	
	If permeate water purity or flow rate is not adequate and does not meet predicted performance levels. Please contact product support for more information: https://www.elgalabwater.com/support/product-support	



9. TROUBLESHOOTING

Issues	Action	
Message Displayed - Sanitization Reminder	1. Complete a Sanitization Procedure. (See section 8.6 / 8.7 Sanitization Procedures)	
Message Displayed - UV Lamp Change Reminder	1. Complete a UV Lamp change. (See section 8.1 Replacing Ultraviolet Lamp)	
	2. Reset UV Lamp reminder. (See section 8.1 Replacing Ultraviolet Lamp)	
Message Displayed - High Temperature Water Alarm	 Adjust the alarm points settings via Main Menu. or 2. Check the feedwater temperature. 3. Press the Dispense button to draw some water through the system, swith OFE and restart. 	
	1. Check the mesh feedwater inlet filter for blockages. (See section 8.3 Cleaning the Inlet Filter) or	
Reduced Dispensing Flow Rate	2. Replace the Point-Of-Use Filter. (See section 8.5 Replacing the Vent Filter LC291 and Point-Of-Use Filter LC134 - LC145 - LC197)	
	or 3. Check the feedwater valves are fully opened. (See section 8.3 Cleaning the Inlet Filter)	
	1. Check the alarm points settings via Main Menu. (See section 10.1 Set Alarm Points)	
Message Displayed - Dispense Water Purity Alarm	or 2. Replace the Purification Pack LC292 and wait for the cycle to complete. (See section 8.2 Replacing the Purification Pack LC292)	
	or 3. Replace the Ultraviolet Lamp LC210. (See section 8.1 Ultraviolet Lamp LC210)	
	(Water Purity Alarm is active during commissioning, see quick start)	



IO. CUSTOMIZE OPERATION

10.1 Set Alarm Points

- 1. Hold and press accept \leftarrow to enter the main menu (Fig.1).
- 2. Scroll $\wedge + \vee$ and select Set Alarm Points, press accept \leftarrow to confirm (Fig.2).
- 3. Scroll $\wedge + \vee$ and select a feature, press accept \leftarrow to confirm (Fig.3).
- 4. Adjust the valve by scrolling $\wedge + \vee$, press accept \leftarrow to save (Fig.4).
- 5. Normal operation will now resume. (Fig.5).



(Fig.1)



(Fig.3)

(Fig.4)

(Fig.5)

10.2 Advanced Data Logging

- 1. Insert a USB stick into the USB port (Fig.1).
- 2. Scroll $\wedge + \vee$ and select 'Start Data Logging', press accept \leftarrow to confirm (Fig.2).
- 3. Scroll \wedge + \vee and select an interval data recording time, press accept \leftarrow to confirm (Fig.3).
- 4. To finish recording press the down \checkmark and remove the USB. (Fig.4).
- 5. Downloaded Data files can be viewed on your PC or Mac.



10.3 Software Update

- 1. Download the latest software file onto a blank USB at: www.elgalabwater.com/customize
- 2. Scroll \wedge + \vee and select 'Software Update', press accept \leftarrow to confirm (Fig.2).
- 3. Progress will be displayed on screen (Fig.3).
- 4. Once completed PURELAB QUEST will power OFF, hold and press ← accept to power ON (Fig.4).
- 5. Software update is now complete, normal operation will resume (Fig.5).





(Fig.1)



(Fig.5)



II. TECHINICAL SPECIFICATIONS

Techinical Specifications

Feedwater Specifications				
Model	PURELAB [®] QUEST RO/DI/UV - PQDIUVM1	PURELAB [®] QUEST RO/DI - PQDIXXM1		
Water Source	(Potable Water Source)	(Potable Water Source)		
Conductivity	<2000µS/cm	<2000µS/cm		
	(high conductivity feedwaters may lower purification pack life and raise Type 3 water conductivity)	(high conductivity feedwaters may lower purification pack life and raise Type 3 water conductivity)		
Hardness	<350 ppm as CaCO 3	<350 ppm as CaCO 3		
Free Chlorine	<0.05 ppm Cl ₂	<0.05 ppm Cl ₂		
Chloramine	<0.02 ppm Cl 2	<0.02 ppm Cl ₂		
Total Chlorine	<0.05 ppm Cl ₂	<0.05 ppm Cl ₂		
Silica	<30 ppm SiO ₂	<30 ppm SiO 2		
Carbon Dioxide (CO ₂)	<30 ppm (recommended <20 ppm)	<30 ppm (recommended <20 ppm)		
Fouling Index	<10	<10		
Iron / Manganese	<0.5 ppm Fe/Mn	<0.5 ppm Fe/Mn		
TOC (Total Organic Carbon)	Recommended <2 ppm	Recommended <2 ppm		
Environment				
Temperature	4 - 40 °C (recommended 10 - 25 °C)	4 - 40 °C (recommended 10 - 25 °C)		
Humidity	Non-condensating. Humidity Max 80%	Non-condensating. Humidity Max 80%		
Surroundings	Clean laboratory environment	Clean laboratory environment		
Feedwater Pressure and Flow Rates				
Feedwater Flowrate	Up to 75 L/hr	Up to 75 L/hr		
Mains Drain Requirements	>90 L/hr	>90 L/hr		
Minimum Inlet Pressure	2 bar (30 psi)	2 bar (30 psi)		
Maximum Inlet Pressure	6 bar (30 psi)	6 bar (30 psi)		
Optimum Inlet Pressure	Regulate to 4 bar (LA512)	Regulate to 4 bar (LA512)		
Pipe Connections				
Inlet	8mm (5/16) OD Tube	8mm (5/16) OD Tube		
Outlet	8mm (5/16) OD Tube	8mm (5/16) OD Tube		
Drain	8mm (5/16) OD Tube	8mm (5/16) OD Tube		
Reservoir Outlets	8mm (5/16) OD Tube	8mm (5/16) OD Tube		
Reservoir Overflow	8mm (5/16) OD Tube	8mm (5/16) OD Tube		



II. TECHINICAL SPECIFICATIONS

Techinical Specifications				
Dimensions and Weights	PURELAB [®] QUEST RO/DI/UV - PQDIUVM1	PURELAB [®] QUEST RO/DI - PQDIXXM1		
Dispatch Weight	422mm / 16.6 In 422mm / 16.6 In F CO C / We IS 232mm / 9.133 In 19.6KG	422mm / 16.6 In 422mm / 16.6 In 5 6 7 7 7 7 7 7 7 7 7 7 7 7 7		
Operational Weight	23KG	21.4KG		
Installation	Worktop or Wall Mounting	Worktop or Wall Mounting		
Electrical Requirements	· · ·			
Mains Input	100-240 VAC, 50-60 Hz	100-240 VAC, 50-60 Hz		
Power Required (Excluding Pump and UV)	24 V DC	24 V DC		
Power Consumption	120 VA	120 VA		
Noise Output	dBA - < 40	dBA - < 40		
Type 1 (Dispenser Tip)				
Resistivity	18.2MΩ.cm@25°C	18.2MΩ.cm@25°C		
Dispenser flow-rate	up to 1.2l/min	up to 1.2l/min		
Available Volume*	Up to 10l/day	Up to 10l/day		
ТОС	<5ppb	<30ppb		
Bacteria (LC134 / LC145 / LC197 POU filter Fitted)	<0.1cfu/ml	<0.1cfu/ml		
Endotoxin (LC197 POU filter fitted)	<0.001Eu/ml	<0.001Eu/ml		
RNases	<1pg/ml	N/A		
DNases	<5pg/ml	N/A		
рН	Effectively Neutral	Effectively Neutral		
Particulates (LC134 / LC145 / LC197 fitted)	0.2µm filtration	0.2µm filtration		
Type 2 Water Outlet (Port: 4)				
Resistivity	>1.0MΩ.cm@25°C	>1.0MΩ.cm@25°C		
TOC (ppbC)	<50ppb	<50ppb		
Bacterial TVC	<100 (cfu/ml)	<100 (cfu/ml)		
Available Volume*	Up to 10L/day	Up to 10L/day		
Type 3 Water Outlet (Port: 5)				
Conductivity (µS/cm)**	<20 (µS/cm)	<20 (µS/cm)		
	<200ppb	<200ppb		
	<1000 (cfu/ml)	<1000 (cfu/ml)		
Darticulator and Pactoria rejection	>>6%	>40%		
Provide the section (MM/s 200 Dr.)	> 2004	> 200/		
Product flow***	>>>>%	>>>>%		
Available volume	υρ το βυι/day	υρ το βυί/day		

* Available Volume of Type 1 and 2 water, combined; Increased use will reduce purification pack life.

**Subject to suitable feed-water purity (see ionic rejection) and system maintenance.



I2. WARRANTY / CONDITIONS OF SALE

ELGA Veolia® is a trading name of VWS (UK) Ltd.

General Limited Warranty

VWS (UK) Ltd. warrants the products manufactured by it against defects in materials and workmanship when used in accordance with applicable instructions for a period of one year from the date of shipment for the products. VWS (UK) Ltd. MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

The warranty provided herein and the data, specifications and descriptions of the VWS (UK) Ltd. products appearing in VWS (UK) Ltd. published catalogues and product literature may not be altered except by express written agreement signed by an officer of VWS (UK) Ltd. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and, if given, should not be relied upon.

In the event of a breach of the foregoing warranty, VWS (UK) Ltd. sole obligation shall be to repair or replace, at its option, any product or part thereof that proves to be defective in materials or workmanship within the warranty period, provided the customer notifies VWS (UK) Ltd. promptly of any such defect. The exclusive remedy provided herein shall not be deemed to have failed of its essential purpose so long as VWS (UK) Ltd. is willing and able to repair or replace any nonconforming VWS (UK) Ltd. product or part. VWS (UK) shall not be liable for consequential, incidental, special or any other indirect damages resulting from economic loss or property damage sustained by any customer from the use of its products.

Water System Limited Warranty

VWS (UK) Ltd. warrants the water systems manufactured by it, BUT EXCLUDING MEMBRANES AND PURIFICATION PACKS, against defects in materials and workmanship when used in accordance with the applicable instructions and within the operating conditions specified for the systems for a period of one year from the earlier of:

the date of installation, or the 120th day following the date of shipment.

VWS (UK) LTD. MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. The warranty provided herein and the data, specifications and descriptions of the VWS (UK) Ltd. systems appearing in VWS (UK) Ltd. published catalogues and product literature may not be altered except by express written agreement signed by an officer VWS (UK) Ltd. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorised and, if given, should not be relied upon.

In the event of a breach of the foregoing warranty, VWS (UK) Ltd. sole obligation shall be to repair or replace, at its option, any product or part thereof that proves to be defective in materials or workmanship within the warranty period, provided the customer notifies VWS (UK) Ltd. promptly of any such defect. The cost of labour for the first ninety (90) days of the above warranty period is included in the warranty; thereafter, labour cost shall be at the customer's expense. The exclusive remedy provided herein shall not be deemed to have failed of its essential purpose so long as VWS (UK) Ltd. is willing and able to repair or replace any nonconforming VWS (UK) Ltd. system or component part. VWS (UK) Ltd. shall not be liable for consequential, incidental, special or any other indirect damages resulting from economic loss or property damage sustained by any customer from the use of its process systems.

Products or components manufactured by companies other than VWS (UK) Ltd. or its affiliates

("Non - VWS (UK) Ltd. products") are covered by the warranty, if any, extended by the Product manufacturer. VWS (UK) Ltd. hereby assigns to the purchaser any such warranty; however VWS (UK) LTD. EXPRESSLY DISCLAIMS ANY WARRANTY WHETHER EXPRESSED OR IMPLIED, THAT THE NON - VWS (UK) Ltd. PRODUCTS ARE MERCHANTABLE OR FIT FOR A PARTICULAR PURPOSE.

NOTICE

VWS (UK) Ltd. is constantly striving to improve its products and services. Consequently, the information in this document is subject to change without notice and should not be construed as a commitment by VWS (UK) Ltd.

Also VWS (UK) Ltd. assumes no responsibility for any errors that may appear in this document. This manual is believed to be complete and accurate at the time of publication. In no event shall VWS (UK) Ltd. be liable for incidental or consequential damages in connection with or arising from the use of this manual.

VWS (UK) Ltd. warrants its products against defects in materials and workmanship as described in the Warranty statement on the preceding pages.







The Labwater Specialists

ELGA Veolia is the global laboratory water brand name of Veolia. The information contained in this document is the property of VWS (UK) Ltd, trading as ELGA Veolia, and is supplied without liability for errors or omissions. © VWS (UK) Ltd. 2020 - All right reserved. ELGA[®], PURELAB[®] are registered trademarks of VWS (UK) Ltd.

