

	EPA SPECIFICATION SHEET	EPA/E/1116-115
		REVISION 01

The following document has now been written and approved. Information relating to the document supply format is stated below. This document is intended to be circulated within ELGA ONLY.

SECTION 1					FRONT COVER				
Part Number		INST41771		Version 01					
Title		PF3-M3 – PURELAB flex 3 - Instruction sheet							
Manual type		Quick Start Manual			<input checked="" type="checkbox"/>	Operator Manual			<input type="checkbox"/>

HEALTH AND SAFETY SHEET




For the Full Operators Manual please scan this QR code or copy this link into your search bar: <https://www.elgalabwater.com/operating-manuals>

ELGA LabWater, Lane End Business Park, Lane End, High Wycombe, HP14 3BY, UK

Tel: +44 (0) 203 567 7300
Fax: +44 (0) 203 567 7305
E-mail: info@elgalabwater.com

For any technical queries please contact techsupport@elgalabwater.com

WARNING! WARNINGS ARE GIVEN WHERE FAILING TO OBSERVE THE INSTRUCTIONS COULD RESULT IN INJURY OR FATALITY!

CAUTION! Cautions are given where failure to observe the instructions could result in damage to the equipment, associated equipment and processes.

WARNING! THINK BEFORE LIFTING! CONSIDER CORRECT LIFTING TECHNIQUES TO AVOID INJURY!

Environmental Requirements

The system should be installed on a flat, level surface, in a clean, dry environment. The system is designed to operate safely under the following conditions:

- Indoor use only
- Altitude up to 2000m
- Temperature Range 5°C - 40°C
- Storage Conditions 2°C - 50°C
- Maximum Relative Humidity 80% @ 31°C decreasing linearly to 50% @ 40°C non-condensing
- The system is in installation Category II, Pollution Degree 2, as per EN 61010-1.
- Noise Levels - <45 dBA

Display Handset

CAUTION! The display handset is not designed to be submerged in water.

Electricity

The appliance coupler (mains lead) or power supply connected to the rear of the unit can be removed to isolate the power supply. If access to this is restricted then it is recommended that access to the supply socket is easily available to disconnect the electrical supply.

WARNING! ONLY USE THE APPLIANCE COUPLER (MAINS LEAD) AND POWER SUPPLY PROVIDED. THE USE OF THESE WILL ENSURE ADEQUATE EARTH PROTECTION IS PROVIDED. IF THE EQUIPMENT IS USED IN A MANNER NOT SPECIFIED BY ELGA, THE PROTECTION PROVIDED BY THE EQUIPMENT MAY BE IMPAIRED. POSITION THE POWER SUPPLY SO THAT IT CANNOT COME INTO CONTACT WITH WATER.

Ultraviolet Light

WARNING! UNDER NO CIRCUMSTANCES SHOULD THE UV LAMP BE CONNECTED AND ACTIVATED WHEN OUTSIDE THE HOUSING. EXPOSING COULD CAUSE SERIOUS INJURY TO EYES AND SKIN. ENSURE THE UV LAMP IS DISPOSED OF IN ACCORDANCE WITH LOCAL REGULATIONS.

PURELAB FLEX INST41771 Ver 1 - 02/24
1

SECTION 2

EXTERNAL CARTON LABELLING

SECTION 2-1

Generic position of carton labels & Carton label details

Carton Product Label



Carton UK EPA
Establishment
Number Label

**EPA Est. 091083-GBR-001
PURELAB**

ELGA **UK** **CE** 

Model PF3-M3



Date of Manufacture 18/01/2024

Serial No. FLC00000003



Maximum Pressure

Inlet 6.0/90 (bar / psi)


Working 6.0/90 (bar / psi)

Electrical 24Vdc-120VA
Supply





REFER TO OPERATING INSTRUCTIONS
SE REFERER AU MANUEL D'OPERATION
IN DER BETRIEBSANLEITUNG NACHSEHEN

ELGA® is the global laboratory water brand of
Veolia Water Solutions & Technologies.
Tel: +44 203 567 7300, Web: www.elgalabwater.com
MADE IN UK

	EPA SPECIFICATION SHEET	EPA/E/1116-115
		REVISION 01



SECTION 3	PRODUCT LABEL DETAIL
------------------	-----------------------------

SECTION 3-1	Generic position of units labels & Unit label details.
--------------------	---



AUTHORISATION - RESEARCH AND DEVELOPMENT DEPARTMENT
--

This document is released on behalf of ELGA and remains the property of VWS (UK) LTD.

Compiled By: 	Date: 17/01/2025
Approved: 	Release Date: 17/1/25