

### The PURELAB<sup>®</sup> Range

### UNINTERRUPTED DISCOVERY

The laboratory water purification solutions for your research needs.



A /

### Contents

About ELGA	05
Your Water	08
Technologies	10
Product Range	12

#### **Product Overview**

Quest	14	Pharma Compliance	28	
Flex 1	18	Chorus 2+	30	
Flex 2	19	Chorus 2	32	
Flex 3	20	Chorus 3	34	
Flex 3+	21	Storage Reservoir	36	
Chorus 1	24	Remote Dispenser	40	
Chorus 1 Complete	26	Halo Dispenser	42	

Product Specification			45	
Quest	46	Chorus 1 Complete	51	
Flex 1 & 2	48	Pharma Compliance	52	
Flex 3 & 3+	49	Chorus 2+	53	
Chorus 1	50	Chorus 2 & 3	54	

#### Index - Product Part Numbers

57





### ELGA. We are the LabWater Specialists

We are the specialists in the engineering, service & support of water purification systems. Our unrivalled product design has achieved international recognition and awards. We have been working with scientists since 1937 to guarantee ultrapure and pure water for their experiments and laboratory work.

### Why choose ELGA as your laboratory water partners?



### Customer-focused – what we create is for our customer

Our commitment to developing and providing you with ultrapure water means that you can focus and concentrate on obtaining accurate and reliable results.



### Innovative – the keystone of our thinking

Our UK R&D facilities are always looking to provide products dedicated to supplying you with the right water quality for your application.



### Sustainable – at the forefront of all of our activities

Our products are designed to have the lowest possible impact on the environment at all stages: manufacture, in service and at end of life. We can calculate the carbon value of all our products through their lifetime.



#### British Engineering and Manufacturing – the standard for all our products

All our systems are manufactured in the UK and we are accredited to ISO:9001 and ISO:14001 standards.

Our technical service teams support our global science & healthcare customers with specialist expertise. Global digital performance monitoring ensures laboratory work is uninterrupted, and a global supply chain supports clients from regional centres on every continent.

#### Our Awards



design effectiveness award winner

reddot design award

winner 2011







### The choice *you* need, for *your* lab

We understand how important it is for you to obtain a choice of water qualities that range from Reverse Osmosis grade for simple routine washing and rinsing, through to ultrapure water for the most critical applications. The PURELAB product range has a wide variety of water purification systems that will meet any one of your requirements for water quality.

## What are the different types of water quality and what are they used for?

#### Type I Water

Often referred to as ultrapure water, this grade is required for some of the most water-critical applications such as HPLC (High-Performance Liquid Chromatography) mobile phase preparation, blanks and sample dilution for other key analytical techniques such as GC (Gas Chromatography), AAS (Atomic Absorption Spectrophotometry), and ICP-MS (Inductively Coupled Plasma Mass Spectrometry) as well as molecular biology applications.

#### Type II Water

Is the grade for general laboratory applications such as media preparation, pH solutions, buffers and used as a feed source to Type I systems.

#### Type III Water

Is the grade of water recommended for non-critical work which may include glassware rinsing, water baths, autoclave, environmental chambers as used to feed Type I systems.

#### There is a PURELAB® product to suit your specific needs

TYPE I*	TYPE II	TYPE III
PURELAB Quest	PURELAB Quest	PURELAB Quest
PURELAB Flex 1**, 2, 3 & 3+	PURELAB Chorus 2 RO/EDI/UV	PURELAB Chorus 3
PURELAB Chorus 1 Analytical Research*	PURELAB Chorus 2 RO/DI/UV	PURELAB Flex 1**
PURELAB Chorus 1 Life Science*	PURELAB Chorus 2 RO/DI	
PURELAB Chorus 1 General Science*	PURELAB Flex 1**, 3 & 3+	
PURELAB Chorus 1 Complete		
DUDELAD Destroy Commissions*		

PURELAB Pharma Compliance\*

\*PureSure Technology/in-line filtration/variety of purification packs . \*\*Dependent on feed water.

# Our comprehensive range of products

















# Technologies

The various technologies used in ELGA equipment are able to remove impurities from water down to extremely low levels; some technologies focus on specific contaminants while others have a broader spectrum of targets. To achieve the correct water purity for a particular application, in a cost effective manner, technologies must be arranged in combination and their operation optimized.

#### **Activated Carbon**

Contains a maze of tiny pores with sizes ranging from <0.1 - 2.9 nm and a surface area of about 1000 square meters per gram. The nature of this surface allows adsorption of organic impurities from the water and catalytic decomposition of free chlorine and, more slowly, chloramines.

It is applied in:

- Pre-treatment cartridges
- Composite vent filters
- Final purification cartridges

#### **Microporous Depth Filters**

Pre-filtration using microporous depth filters, provides an entrapment/adsorption barrier for the removal of large suspended particles and some colloids from the water entering the ELGA purification process. Typically rated at  $5 - 10 \mu m$ and combined with an activated carbon treatment, these filters act to protect subsequent RO systems from fouling and blockage.

#### Small and Large Organic Molecules



#### **Reverse Osmosis**

Reverse osmosis (RO) is a process where pressure is used to push water through a membrane filter in a cross-flow fashion. RO membranes are extremely fine filters that reject water contaminants less than 1nm in diamenter. Typically >95% of ionic impurity, most organic impurity, and nearly all particulates, bacteria and bio-molecules are removed from the permeate water; these are carried out of the RO module in a waste or concentrate water stream.



## Technologies

#### Ion Exchange

lon exchange resins are often used as part of a final treatment step. Single-use purification packs typically use a mixture of ion-exchange resins and other media. When used to deionize water, charged impurities are retained on these resins, while H+ and OH- ions (which combine to form water) are released to replace them. This is a highly efficient process and can remove ions to give product-water resistivity of up to 18.2 MΩ.cm.

#### Electrodeionization

Electrodeionization (EDI) combines ion-exchange resins and ion-selective membranes, which are used to move ionic impurities into a waste or concentrate water stream leaving purified productwater. As impurities leave via the concentrate water stream, their build-up does not exhaust the resin, and therefore prolongs resin lifespan. A single EDI unit may operate for many years before a replacement is required. Typically product water resistivity of >10 M $\Omega$ .cm is consistently achieved using this process. This technology can be used as an alternative to single-use purification cartridges.

#### Ultraviolet (UV) Light

Treatment of water with UV-C light is used to photo-oxidise organic impurities and/or inactivate micro-organisms. Photo-oxidation of organic impurities results in polar or charged species that can subsequently be removed by ion-exchange processes. Typically the UV lamp forms part of a 'polishing' treatment loop including ion-exchange, through which water is repeatedly circulated to maintain quality. Water with Total Organic Carbon (TOC) of <5 ppbC and bacteria at <1 CFU/ml can be achieved in ELGA products that use this approach.

#### **Sub-micron Filtration**

Sub-micron filtration, including micro, ultra-micro and ultra filters (30–3000 nm) are used as part of a 'polishing' loop or at the point-of-use. Fine filtration is applied to remove bacteria (live or dead) and biologically active molecules. These absolute filters have pores smaller than their intended target and can retain the impurity while allowing water to pass through. Impurities that are removed by sub-micron filtration, include bacteria, colloids, enzymes, endotoxins and particulates.

#### **Environmental improvements**

New Aquaporin membrane technology provides greater recovery with reduced membrane surface area, reducing both water and component waste.



## **PURELAB**<sup>®</sup>

**Product Range** 



### PURELAB<sup>®</sup> Quest

The only purifier on the market that dispenses all 3 types of science ready water from a compact, economical and easy to use system. The Quest water solution system provides laboratory water directly from a tap water input.

#### AT A GLANCE

#### **Overview:**

- Key Features:
- Type I 18.2 MΩ.cm (Up to 10 Litres per day)
- Type II > 1 MΩ.cm (Up to 10 Litres per day)
- Type III 20 μS.cm (Up to 30 Litres per day)
- Compact
- Connected
- Sustainable
- Cost effective
- Proven reliability which is tested, proven and unrivalled
- Effortless, ingenious, intuitive

#### **Ideal Applications:**

- Type I Water HPLC GC-MS AA/ICP-OES IC
- Type II Water Preparing and diluting buffers and reagents
- Type III Water Glassware rinsing Water baths

Molecular biology e.g. DNA Sequencing and PCR

Tissue culture media pH solutions

Autoclave feeds

#### Compact

232 mm wide. 510 mm high. Wall mountable; saving valuable lab space.

#### **Sustainable**

Made from more than 85% reclaimed materials.\* Designed with long lasting consumables.

#### **Cost-Effective**

3 types of water from one system. Space saving design means a more efficient lab and team.



\*Excludes materials in contact with the purification process.

*3 types of water from one system. Space saving design means a more efficient lab and team.* 

UREL

ELG/



Troum

### The PURELAB® flex Range









flex 2



flex 3+



### PURELAB® flex 1

#### **Simplicity and Elegance**

A versatile dispensing and monitoring system when connected to a reservoir or distribution loop, while also functioning seamlessly as a straightforward deionization system.

#### AT A GLANCE

#### **Overview:**

- Type I
- Type II
- Type III

Dependant on use of DI pack Litres per day: Dependent on feed water

#### **Key Features:**

- Flexible dispensing
- Customise settings
- Full re-circulating
- Fully Integrated filtration

#### **Ideal Applications:**

- General Lab
- Type II Applications
- Dispensing with DI Pack

#### **Customized Settings**

Be in control of your PURELAB flex by customizing the settings to suit your application.

#### Simplicity

Simple to install, operate and maintain with prioritized information displayed at all times (system status, alarm) for absolute confidence as you dispense.

#### Easy to Use

Ergonomic design with features including auto volumetric dispense and height adjustability.

#### Data Capture

Data capture via USB for system performance validation and software updates.

#### Space Saving

Space saving and compact dispenser which can be placed on the bench or wall mounted.

#### Sustainability

Manufactured using up to 70% reclaimed materials.

See specification on page 48

### PURELAB® flex 2

#### Designed for the laboratory of today

Ensures precise, adaptable, and user-friendly performance. This acclaimed system generates ultrapure type I (18.2  $\Omega$ .cm) water from a pre-purified source, making it perfect for analytical and life science purposes.

#### AT A GLANCE

#### **Overview:**

- **Type I** 18.2 MΩ.cm
- Litres per day: Dependent on feed water

#### **Key Features:**

- Real-time TOC
- Fully re-circulating
- Customize settings
- Adjustable dispensing

#### **Ideal Applications:**

- Mass Spectrometry
- Molecular Biology
- Electrochemistry
- Atomic Spectroscopy
- Liquid Chromatography
- Gas Chromatography
- Immunochemistry
- Spectrophotometry
- Media / Buffer Prep
- General Chemistry

#### **Full Spectrum UV Treatment**

Full recirculation through the UV lamp and purification pack right to the point of use for peace of mind.

#### **Intuitive Flexible Dispense**

Clear water purity display for absolute confidence as you dispense.

#### **Real-time TOC Monitoring**

Provides complete confidence in organic purity by reducing the level of organics for critical applications.

#### Sustainability

Utilising PURELAB flex new ECO MODE reduced power consumption by 17% and manufactured from up to 70% reclaimed material.

#### **Data Capture**

Data capture via USB for system performance validation and software updates.

#### Space Saving

Space saving and compact dispenser which can be placed on the bench or wall mounted.

### PURELAB® flex 3

#### **Power and flexibility**

A small unit with big powerful capabilities. The PURELAB flex 3 is the ultimate system providing Type 1 ultrapure water from potable tap water in one single unit.

#### AT A GLANCE

#### **Overview:**

- **Type I** 18.2 MΩ.cm
- Litres per day: < 10

#### Key Features:

- Real-time TOC
- Fully re-circulating
- Adjustable dispensing

#### **Ideal Applications:**

- Mass Spectrometry
- Molecular Biology
- Electrochemistry
- Atomic Spectroscopy
- Liquid Chromatography
- Cell Culture

- Gas Chromatography
- Immunochemistry
- Spectrophotometry
- Media / Buffer Prep
- General Chemistry

#### Space Saving Design

The compact unit can be placed on the bench or wall mounted and has an integrated 7 Litre reservoir filled by a 20 l/hr RO membrane, ensuring that water is always available.

#### **Full Recirculating**

Ensuring the highest microbial purity and guaranteeing pure water, as recirculation of the water occurs from the reservoir right to the point-of-use.

#### **Flexible Dispenser**

The intuitive dispenser offers a clear display of the water purity for absolute confidence as you dispense.

#### **Real-time TOC Monitoring**

Provides complete confidence in organic purity and clear display at all times. The final quality sensor is placed at the entry of the flexible dispenser giving you peace of mind.

#### Sustainability

Manufactured using up to 55% reclaimed materials; improvements in energy efficiency of up to 55%, and reduced water waste up to 20%<sup>\*</sup>.

\*Compared to previous generation PURELAB Flex 3.

#### Data Capture

Data capture via USB for system performance validation and software updates.

See specification on page 49

### **PURELAB**<sup>®</sup> flex 3+

#### **Flexible elegance**

A unique automated solution providing Type I (ultrapure water) from potable water. Couples directly with analytical chemistry systems delivering complete automation - all in one single unit.

#### AT A GLANCE

#### **Overview:**

#### **Key Features:**

- Litres per day: <10
- Type I 18.2 MΩ.cm Connects directly to multiple analytical chemistry systems from multiple brands
  - On-demand continuous ultrapure water to analysers
  - Full automation of water production & delivery
  - Real-time TOC monitoring
  - Fully re-circulating
  - Adjustable dispensing

#### **Ideal Applications:**

- Analytical Chemistry systems
- Mass Spectrometry
- Molecular Biology
- Electrochemistry
- Atomic Spectroscopy
- Liquid Chromatography
- Cell Culture
- Gas Chromatography
- Immunochemistry
- Spectrophotometry

#### **On-demand Water Supply**

Provides consistent & continuous supply of ultrapure water to analysers to reduce failed runs and loss of samples, reagents and analyser components/accessories.

#### **Space Saving Design**

The compact unit can be placed on the bench or can be wall mounted and has an integrated 7 Litre reservoir filled by a 20 l/hr RO membrane, ensuring that water is always available.

#### **Full Recirculating**

Ensuring the highest microbial purity and guaranteeing pure water, as recirculation of the water occurs from the reservoir right to the point-of-use.

#### **Real-time TOC Monitoring**

Provides complete confidence in organic purity and clear display at all times. The final quality sensor is placed at the entry of the flexible dispenser giving you peace of mind.

#### Data Capture

Monitor water purity via USB datalog ensuring analytical performance during test batch.

#### **Sustainability**

Manufactured using up to 52% reclaimed materials; improvements in energy efficiency of up to 55%, and reduction of water to drain by up to 20%\*.

\*Compared to previous generation PURELAB Flex 5.

### The Chorus Range

**Chorus 1** 





Pharma Compliance



Chorus 2+



Chorus 2

Chorus 3







### PURELAB<sup>®</sup> Chorus 1

#### Flexible. Configurable. Simple.

Delivering the ultimate in water purity for absolute confidence in your results. When you require the ultimate in water purity, PURELAB Chorus 1 provides the perfect solution. Consistently delivering water purity of 18.2 M $\Omega$ .cm (Type I+/I) and underpinned by the advanced PureSure<sup>®</sup> deionization system.

#### AT A GLANCE

#### **Overview:**

- Key Features:
- Type I+ Water
  18.2 MΩ.cm
- Litres per day: Depends on feedwater
- Real time TOC
- Fully re-circulating
- Integrated filtration
- Multiple dispensing
- PureSure technology

#### **Ideal Applications:**

- Mass Spectrometry
- Ion Chromatography
- Ultra trace Analyses
- Qualitative Analyses
- Gas Chromatography
- Molecular Biology
- Cell Cultures
- Electrochemistry
- Immunochemistry
- Atomic Spectroscopy

#### **Advanced PureSure Deionization**

A twin-bed ion-exchange process with inter-stage resistivity monitoring allows retention of any impurity released during exhaustion of a primary cartridge, by the secondary polishing cartridge. This method gives guaranteed, optimum, product-water quality, advanced warning of consumable change, and extended consumable service life.

#### **Fully Recirculating**

Ensuring microbial purity and guaranteeing pure water at the point-of-use.

#### **Real-time TOC Monitoring**

Provides complete confidence in organic purity.

#### **Integrated Filtration**

Ultrafiltration or microfiltration filters out endotoxins, proteins, nucleases and particulates. 185nm wavelength breaks organic compounds and 254nm wavelength sterilises bacteria and viruses.

#### Data Capture

Data capture via USB for system performance validation and software updates.



Chorus 1 enables you to focus on attaining accurate results while ensuring an uninterrupted work flow

PURELAB

ELGA

Ċ



### PURELAB<sup>®</sup> Chorus 1 Complete

#### One complete solution for the laboratory.

PURELAB Chorus 1 Complete provides a complete solution from potable tap water supply to ultrapure water. With its easy to use ergonomic design, new and improved parameter readings and dispensed with confidence directly from the system or from a choice of additional Halo Dispensers.

#### AT A GLANCE

#### **Overview:**

- Type I Water
  18.2 MΩ.cm
- Litres per day: Up to 480
- **Key Features:**
- Tap-to-ultrapure
- Real time TOC\*
- Fully recirculating
- Multiple dispensing

#### **Ideal Applications:**

- Mass Spectrometry
- Molecular Biology
- Electrochemistry
- Atomic Spectroscopy
- Liquid Chromatography
- Gas Chromatography
- Immunochemistry
- General Laboratory
- Spectrophotometry

#### Tap to Type 1 water

Generate laboratory-grade Type 1 ultrapure water directly from your tap with our all-in-one purification system.

#### **Realtime TOC monitoring\***

Provides complete confidence in organic purity, with optional Inhibit mode for ensured water quality confidence.

#### **Fully Recirculating**

Recirculation of purified water through our modular reservoir to maintain consistent peak water purity at  $18.2 \text{ M}\Omega$ .cm.

#### Aquaporin<sup>™</sup> membrane

ELGA's incorporation of new Aquaporin membrane technology provides 20 L/hr with reduced membrane surface area, reducing component waste.

#### Type I and type II water

Through effective recirculation technology, the Chorus 1 Complete simultaneously produces Type 1 ultrapure water and stores Type 2 purified water within the reservoir, offering dual-grade water production in a single, efficient system.

#### Data Capture

Data capture via USB for system performance validation and software updates.

\*Dependent on model.



See specification on page 51

### **PURELAB**<sup>®</sup> **Pharma Compliance**

PURELAB Pharma Compliance is packed with innovative features optimised for QC laboratories. It offers all the necessary purification, software capabilities, qualification documentation and online support services necessary for the validation required to meet Good Manufacturing Practice (GMP). The ultrapure water system meets your laboratories needs, with minimal intervention, improved efficiencies and a contemporary style to compliment any modern laboratory.

The Pharma Compliance includes secure software providing the evidence to meet FDA and EU standards for digital record keeping. Purification processes fulfil United States Pharmacopeia standards 643 and 645. Qualification documents provided enable the unit to support labs subject to validation under GxP.

#### AT A GLANCE

#### **Overview:**

• Type I Water - 18.2 MΩ.cm

#### **Key Features:**

- Smart, intuitive software
- Digital record keeping
- Quality management system
- Meets GMP requirements
- Complies with the TOC verification and water conductivity measurements required under US pharmacopeia 643 and 645.

#### **Ideal Applications:**

- QC labs wanting to meet high level GMP regulations
- QC labs with FDA CFR 21 Cell cultures Part 11 requirements
- Mass spectrometry
- Ion and Gas chromatography
- Ultra trace and qualitative analyses
- Molecular Biology
- Electrochemistry
- Immunochemistry

#### **Digital Record Keeping**

The PURELAB Pharma Compliance's admin functions, security and data integrity capabilities, password accessibility, audit trail procedures & permissions all meet national regulators best practice for data management in QC labs following GMP standards.

#### **Designed to Comply**

Designed to meet FDA, United States Pharmacopeia, EUDRALEX, European Pharmacopoeia and all GMP requirements for Quality Control laboratories.

#### **Cost Effective**

The PURELAB Pharma Compliance delivers the water to support the QC tests necessary to validate drug purity in pharma manufacturing. Consumables are designed to minimise cost and waste.



See specification on page 52

### **PURELAB**<sup>®</sup> Chorus 2+

#### Flexible. Configurable. Simple.

PURELAB Chorus 2+ (RO/EDI/UV) features our patented recirculating EDI technology: the only EDI system on the market that is able to fully recirculate to achieve 10- >15 M $\Omega$ .cm. The PURELAB Chorus 2+ provides additional bacteria and inorganic quality for sensitive analytical and life science applications above that of basic laboratory work. With its simple design and ease of use, water can be measured and dispensed with confidence from the system or from a choice of additional Halo or **Remote Dispenser.** 

#### AT A GLANCE

#### **Overview:**

#### **Key Features:**

- Type II+ Water - 10 MΩ.cm
- Litres per day: upto 480L/day (RO/DI/UV) upto 220L/day (RO/EDI/UV)
- Tap to Type II
- Fully Re-circulating
- Multiple Dispensing

#### **Ideal Applications:**

- Electrochemistry
- Cell Culture
- Spectrophotometry
- Feed to Ultrapure Water

#### **Fully Recirculating**

ELGA's patented fully recirculated EDI provides a constant supply of high purity that guarantees a minimum of 10 M $\Omega$ .cm water at all times.

#### **Ideal for High Volume Labs**

A cost-effective solution for laboratories requiring higher output volumes thanks to the incorporated EDI technology.

#### **Single System Solution**

Perfect single system solution for analytical and life science applications.

#### **Space Saving Design**

Designed to be modular and stackable to save space, whether wall-mounted\* or under the bench. \*If fitted with Halo dispenser

WWW.ELGALABWATER.COM

#### **Reduced Maintenance Times**

Quick and easy replacement of consumables as well as semi-automated sanitization to reduce maintenance times.

• Media / Buffer

Preparation

General Chemistry

#### **Economical**

Optional CO, removal from the purified water (post RO) increasing the life of downstream consumables.

#### **Data Capture**

Data capture via USB for system performance validation and software updates.



See specification on page 53

### PURELAB<sup>®</sup> Chorus 2

#### Modular. Flexible. Reliable.

Reliable delivery of Type II water purity. When Type II water is all you need, then PURELAB<sup>®</sup> Chorus 2 (RO/DI) is the reliable solution with the flexibility to suit your requirements. It dispenses up to 480 Litres of pure water per day from a potable water supply for general laboratory applications.

#### AT A GLANCE

#### **Overview:**

- Key Features:Easy Configurability
- Type II 10 MΩ.cm
- Litres per day: Up to 480
- Modular

#### **Ideal Applications:**

- Stills Replacement
- Buffer Preparation

Washing/Rinsing

- pH solution Preparation Steam Generators
  - Sterilizer Feed

Hydroponics

- Autoclaves
- Sterilizer Feed

General Chemistry

Feed to Type I Polishers

#### Deionization

The Purification Pack contains optimized resin mixes to maximize consumables capacity.

#### Simplicity

Simple to install, operate and maintain with prioritized information displayed at all times (system status, alarm) for absolute confidence as you dispense.

#### **Economical**

Optional CO<sub>2</sub> removal from the purified water (post RO) increasing the life of downstream consumables.

Option to reduce water consumption for low hardness feed waters.

#### Modular

Multiple PURELAB Chorus 2 units can feed into one reservoir and systems can be expanded post-installation. As such, the cost of future upgrades is minimized. Duplex systems also guarantee maximum uptime.

#### Data Capture

Data capture via USB for system performance validation and software updates.

#### **Reduced Maintenance Times**

Quick and easy replacement of consumables as well as simple sanitization to reduce maintenance time.

PURELAB Chorus 2 (RO/DI) is the reliable solution with the flexibility to suit your requirements



See specification on page 54

### PURELAB<sup>®</sup> Chorus 3

#### Modular. Flexible. Reliable.

Reliable delivery of Type III water purity. When general laboratory grade water is all you need, then PURELAB Chorus 3 is the reliable solution with the flexibility to suit your requirements. It can also be used as a feed to other ELGA water systems.

#### AT A GLANCE

#### **Overview:**

- Type III Water
  10 MΩ.cm
- Litres per day: Up to 720
- RO Water

**Key Features:** 

- Easy Configurability
- Auto Rinse
- Modular

#### **Ideal Applications:**

- Buffer Preparation
- Washing/Rinsing
- Autoclaves
- General Chemistry
- Hydroponics
- Steam Generators
- Sterilizer Feed
- Feed to Type I polishers

#### Configuration

Ability to configure multiple systems to increase flow rate.

#### Simplicity

Simple to install, operate and maintain with prioritized information displayed at all times (system status, alarm) for absolute confidence as you dispense.

#### **Economical**

Optional CO<sub>2</sub> removal from the purified water (post RO) increasing the life of downstream consumables.

Option to reduce water consumption for low hardness feed waters.

#### **Auto Rinse**

Maintains purity of water during periods of low use.

#### Data Capture

Data capture via USB for system performance validation and software updates.

#### **Reduced Maintenance Times**

Quick and easy replacement of consumables as well as semi-automated sanitization to reduce maintenance time.

#### Modular

Multiple PURELAB Chorus 3 units can feed into one reservoir and systems can be expanded postinstallation. As such, the cost of future upgrades is minimized. Duplex systems also guarantee maximum uptime.



See specification on page 54

### Storage Reservoirs

Our unique range of storage solutions are designed to maintain optimum purity of stored water and provide effective protection against airborne contaminants. They are designed to accommodate PURELAB Chorus water purification systems by maximizing the space in a single integral, compact unit or to sit independently to suit the layout of your laboratory.

#### **Multiple Positioning**

Multiple positioning/mounting options to suit your laboratory layout.

#### **Polyethylene Construction**

Inert opaque polyethylene construction with smooth inner surface.

#### **Dispense Tap**

Positioned to minimize accidental operation or damage (choice of positions).

#### **Advanced Vent Filtration**

Prevents the ingress of airborne bacteria, particulates, organic vapours and CO<sub>2</sub>.

#### **Auto Fill**

Monitoring of reservoir water levels with automated refill ensures purified water is always available.

#### **Hygienic Overflow**

Hygienic overflow in the unlikely event of water system malfunction.

#### **Easy Display**

Direct display of stored purified water on the front of the reservoir for easy identification.



Capacity: 15 Litres Dimensions (mm): 485 (h) x 376 (w) x 347 (d)\* Flow Rate: 6 I/min

\* Please allow for ±8mm tolerance



Capacity: 30 Litres Dimensions (mm): 676 (h) x 376 (w) x 347 (d)\* Flow Rate: 8 I/min



Capacity: 60 Litres Dimensions (mm): 590 (h) x 532 (w) x 534 (d)\* Flow Rate: 10 I/min



Capacity: 100 Litres Dimensions (mm): 805 (h) x 532 (w) x 524 (d)\* Flow Rate: 10 l/min





### **PURELAB**<sup>®</sup> Dispensing Solutions

18.2

A

### PURELAB® Dispenser

A remote dispenser with an in-built water quality monitor – providing maximum reassurance that you are not losing purity as you move to a free-standing solution. Together with PURELAB Chorus & Quest you can quickly and easily create a bespoke water purification solution for your lab. The PURELAB Dispenser enables uninterrupted discovery.

#### AT A GLANCE

#### **Overview:**

#### **Key Features:**

- Type I Water
  18.2 MΩ.cm
- Unrivaled purity in designUnmatched flexibility
- Litres per day: Up to 480
- Effortless use
- Simple & easy installation

#### **Ideal Applications:**

- Life Science applications (cell culture, PCR, genomics, dissolution testing, western blotting)
- Chromatographic techniques
- (HPLC, UHPLC, LC-MS, IC, ICP-MS, AAS )

#### Accuracy

- Flexible dispensing with unrivalled built in purity monitoring sensor
- Delivering proof of water purity
- Repeated volume dispensing
- Improving lab efficiency
- Reproducible results

#### **Streamlined footprint**

- Maximized lab space
- Optimization of storage space
- Up to 4 dispensers per water purifier\*
- \*Dependent on water purifier purchased.

#### Simple use and maintenance

- Plug and Play installation with PURELAB range
- Height adjustable, 180° rotating arm and flexible dispenser arm
- IPX7 rating dispenser handset with menu navigation
- Optional connection to point of use filtration

#### Efficient

- From precise drop by drop up to 2 litres per minute allowing for the right speed including calibration work
- Variable flow rate for easier filling of different sizes of labware
- Volumetric dispensing freeing you to continue working while water is being dispensed
- Profile & locked dispensing allowing you to choose your desired amount of water and prevent accidental flooding



### Halo Dispense Solutions

The PURELAB range offers a variety of dispensing and monitoring solutions to customers enabling the ultimate flexibility when using PURELAB Chorus 1, Chorus 1 Complete, & Quest.

#### **Clear Display**

Prioritized information displayed at all times (system status, alarm) for absolute confidence as you dispense.

#### **Multiple Positioning**

Position the dispenser on the benchtop or independent from the water purification system to optimise your valuable lab space.

#### **Flexible Dispensing**

- Variable flow
- Auto-volume dispense
- Hands free
- Locked dispense
- · Hand-held dispensing
- Profile dispense

#### Halo Glow\*

The unique glow changes colour and flashes alerting you to changes in the system performance.

#### Water monitoring

Water is monitored right to the point of dispense for complete peace of mind.





- Variable flow rate dispense
- Drop-by-drop control
- Locked dispense
  - Purity monitoring to point-of-use
  - Auto volume dispense
  - Profile dispense
- Flexible handset



\*Only for PURELAB Chorus range.

A choice of four dispensing solutions are available to suit different applications, budget and configuration

2

ELG/

ഗ



# **PURELAB**®

Specification



#### Process Flow PURELAB Quest



Specifications	PURELAB Quest UV	PURELAB Quest		
Ultrapure (Type I) water specifications (from fixed dispense tip)				
Resistivity	18.2 MΩ.cm @ 25°C	18.2 MΩ.cm @ 25°C		
Dispense flow rate	Up to 1.2 l/min (0.27 gal)	Up to 1.2 l/min (0.27 gal)		
тос	<5 ppb	<30 ppb		
Bacterial TVC	<0.1 CFU/ml *1	<0.1 CFU/ml *1		
Endotoxin	<0.001 EU/ml *2	<0.001 EU/ml *2		
RNases	<1 pg/ml	n/a		
DNases	<5 pg/ml	n/a		
рН	Effectively Neutral	Effectively Neutral		
Particulates	0.2 $\mu m$ filtration *1	0.2 µm filtration *1		
Recommended daily volume	Up to 10 l/day *3	Up to 10 l/day *3		
Pure (Type II)	) water specifications (Water Out	let port 4)		
Resistivity	>1 MΩ.cm @ 25 °C	>1 MΩ.cm @ 25 °C		
TOC	< 50 ppb	< 50 ppb		
Bacterial TVC	<100 CFU/ml	<100 CFU/ml		
Recommended daily volume	Up to 10 l/day *³ (2.2 gal)	Up to 10 l/day *³ (2.2 gal)		

PURELAB Quest	PURELAB Quest UV		
Operational weight	Operational weight		
21.4 kg (47.2 lbs)	23 kg (50.7 lbs)		
Dime	ensions		
Height	511mm		
Width	232mm		
Depth	421mm		
Insta	Ilation		
Benchtop of	wall-mounted		
Electrical Requirements			
Main Input	100 - 240 VAC, 50 - 60 Hz		
Power Required (Excluding Pump and UV)	24 V DC		
Power Consumption	120 VA		
Noise Output	dBA - <40		
Pipe Co	nnections		
Inlet	8mm (5/16) OD Tube		
Outlet	8mm (5/16) OD Tube		
Drain	8mm (5/16) OD Tube		
Reservoir Outlets	8mm (5/16) OD Tube		
Reservoir Overflow	8mm (5/16) OD Tube		
Envir	onment		
Temperature	4 - 40 °C (recommended 10 - 25 °C)		



#### Process Flow PURELAB Quest UV



RO-permeate (Type III) Water Specifications (Water Outlet port 5)			
Conductivity	<20 µS/cm *4	<20 µS/cm *4	
TOC	<200 ppb *4	<200 ppb *4	
Bacterial TVC	<1000 CFU/ml *4	< 1000 CFU/ml *4	
lonic rejection	>96% *5	>96% *5	
Particulates and Bacteria rejection	>99%	>99%	
Organic rejection (MW > 200 Da)	>99%	>99%	
Production flow	10 l/hr *5	10 l/hr *5	
Recommended daily volume	Up to 30 l/day	Up to 30 l/day	

\*1 When using point of use filters (LC134/LC197)

\*2 When using point of use filter (LC197)

\*3 Available volume of Type I and II water combined; increased use will reduce purification pack life

\*4 Subject to suitable feed water purity (see ionic rejection) and system maintenance

\*5 With feed water pressure at > 4 bar and temperature at 15  $^\circ \! C$ 

Feedwater Specifications			
Water Source	(Portable Water Source)		
Conductivity	< 2000 µS/cm (High conductivity feedwaters may lower purification pack life and raise Type III water conductivity)		
Hardness	<350 ppm as CaCO <sub>3</sub>		
Free Chlorine	<0.05 ppm Cl 2		
Chloramine	<0.02 ppm Cl 2		
Total Chlorine	<0.05 ppm Cl 2		
Silica	<30 ppm SiO <sub>2</sub>		
Carbon Dioxide (CO2)	<30 ppm (recommended < 20 ppm)		
Fouling Index	<10		
Iron/Manganese	<0.5 ppm Fe/Mn		
TOC (Total Organic Carbon)	Recommended <2ppm		

### flex 1 & 2



#### Process Flow PURELAB flex 1







Application	PURELAB Flex 1	PURELAB Flex 1 (with purification pack)	PURELAB flex 2		
Daily volume	>10 Litres	>10 Litres	>10 Litres		
Dispense Flowrate	Up to 2.0 l/min	Up to 2.0 l/min	Up to 2.0 l/min		
Inorganics @25°C	As per feedwater	18.2 MΩ.cm	18.2 MΩ.cm		
Total organic carbon (TOC)	Dependent on	feedwater	<5 ppb		
Bacteria	<0.1 CFU/ml°	<0.1 CFU/ml°	<0.001 CFU/ml <sup>01</sup>		
Bacterial Endotoxin	N/A	<0.001 EU/ml*	<0.001 EU/ml*		
RNase	N/A	N/A	<1 pg/ml		
DNase	N/A N/A <5 p		<5 pg/ml		
$^{\circ}$ With C134/145/197 POU filter/Biofilter $^{*}$ With LC197 Biofilter $^{\scriptscriptstyle 1}$	filter/Biofilter <sup>†</sup> With LC197 Biofilter <sup>1</sup> Dependent on feedwater				
Source	Originally from potable supply, then pretreated. Preferably reverse osmosis (RO) or filtered service deionization (SDI) or distilled.				
Fouling index (max)	<1 for all models				
Free Chlorine	<0.05 ppm max				
тос	N/A				
Carbon dioxide	<0.1 ppm				

Source	Originally from potable supply, then pretreated. Preferably reverse osmosis (RO) or filtered service deionization (SDI) or distilled.		
Fouling index (max)	<1 for all models		
Free Chlorine		<0.05 ppm max	
тос		N/A	
Carbon dioxide		<0.1 ppm	
Silica		<2 ppm	
Particulates	5-10 μm		
Temperature	4-40°C (Recommend 10-15°C)		
Flowrate (maximum requirement)	>2 1/min (0.5 USG)		
Drain requirements	None required		
Feedwater pressure	1.5 bar (22 psi) maximum; Flooded suction minimum		
* Fit LA652 Pressure Regulator where feedwater pressure exceed	A652 Pressure Regulator where feedwater pressure exceeds specified limits		
Dimensions	Height 900-1020mm, Width 236mm, Depth 374mm		m
Weight	10 kg (22 lbs)	10.5 kg (23.1 lbs)	11 kg (24.2 lbs)
Installation	Bench / wall		

## flex 3 & 3+



#### Process Flow PURELAB flex 3



	25217505
Application	PURELAB Flex 3
Daily volume	<10 Litres
Dispense Flowrate	Up to 2.0 l/min
Reverse osmosis make up flow rate at 15°C	20 l/hour
Inorganics @25°C	18.2 MΩ.cm
Total organic carbon (TOC)	<5 ppb
Bacteria	<0.001 CFU/ml°
Bacterial Endotoxin	<0.001 EU/ml*
RNase	<1 pg/ml
DNase	<5 pg/ml
* Dependant on feed water ° With LC134/145/197 PO	U filter/Biofilter <sup>‡</sup> With LC197 Biofilter
Source	Potable tap water
Fouling index (max)	<10
Free Chlorine	<0.5 ppm max
тос	<2 ppm
Carbon dioxide	<30 ppm (recommended <20 ppm)
Silica (recommended max)	<30 ppm
Particulates	-
Temperature	4-40°C (Recommend 10-15°C)
Flowrate (maximum requirement)	Up to 75 l/hr (20 USG)
Drain requirements	<90 l/hr (23 USG)
Feedwater pressure*	6 bar (90 psi) max; 2 bar (30 psi) min
* Fit LA652 Pressure Regulator where feedwater	pressure exceeds specified limits
Dimensions	Height 900-1020mm, Width 236mm, Depth 470mm
Weight	23 kg (57.3 lbs)
Installation	Bench / wall

#### Process Flow PURELAB flex 3+



	Application	PURELAB Flex 3+	
	Daily volume	<10 Litres	
	Delivery flow rate - maximum	Up to 2.0 l/min	
כ	Recirculation flow rate	1 l/min	
	Reverse osmosis make up flow rate	20 l/hour	
	Inorganics (resistivity @25°C)	18.2 MΩ.cm	
2	Organics (TOC) - typical	<5 ppb*	
- 1 1	Direct from internal reservoir	Type III / RO Water	
	Bacteria - typical	<0.001 CFU/ml°	
	Endotoxin	<0.001 EU/ml*	
	DNase	<5 pg/ml	
	RNase	<1 pg/ml	
	* Dependant on feed water ° With LC134,145/197 POU filter/biofilter <sup>1</sup> With LC197 Biofilter		
	Source	Potable tap water	
	Conductivity	<2000 µS/cm²	
_	Contaminant		
	Hardness	<350 ppm as $CaCO_3$	
	Hardness	<0.5 ppm Cl2	
	Silica	<30 ppm SiO2	
	Carbon Dioxide Maximum	<30 ppm (recommended <20 ppm)	
	Fouling index	<10	
2	Iron / Manganese	<0.05 ppm	
	Organics (TOC)	<2 ppm	
	Particulates	N/A	
	Temperature	4 - 40°C (Recommended 10 - 25°C)	
	Flowrate (requirement at 15°C)	Up to 75 l/hr	
	Drain requirements (gravity fall with air gap)	>90L/hr	
-	<sup>2</sup> Purification pack life may vary with feedwaters >1	400 μS/cm	



Process Flow PURELAB Chorus 1 Life Science

- Ultrafiltration particulate filter removal of endotoxin and large molecules for Life Science Applications



Process Flow PURELAB Chorus 1 Analytical Research - Ultra microfiltration to remove smaller particulates for Analytical Research Applications



Application	Life Science Analytical Research General S							
Dispense Flowrate	Up to 2.0 l/min <sup>+</sup>	Up to 2.0 l/min <sup>+</sup>						
Inorganics @25°C	18.2 MΩ.cm	18.2 MΩ.cm	18.2 MΩ.cm					
Total organic carbon (TOC)	1-3 ppb*	3-10 ppb*						
Bacteria	<0.001 CFU/ml°	<0.001 CFU/ml° <0.001 CFU/ml°						
Bacterial Endotoxin	<0.001 EU/ml	<0.001 EU/ml*	<0.001 EU/ml*					
рН	Effectively neutral	Effectively neutral	Effectively neutral					
Particles (filtration)	<0.01 µm	<0.05 µm	<b>0.2</b> μm°					
RNase	<1 pg/ml	<1 pg/ml						
DNase	<5 pg/ml							
Purification pack capacity	Litres to 18.	2 MΩ.cm = 94,100/(μS/cm + (2.3 x ppm CO <sub>2</sub>	))					
* Dependant on feed water – recommended feed <50 ppb TOC ° With LC134/145/197 POU filter/Biofilter 'When connected to Halo, Advanced or flexible 'With LC197 Biofilter								
Source – originally from potable supply, then pretreated	Preferably RO produced by PURELAB Chorus 3 or filtered service deionization (SDI) or distilled. Note: mixed bed or twin bed deionized supplies should be cation limited at exhaustion							
Fouling index (max)	1 for all models. A 5-10 micron membrane prefilter is recommended for all non-RO feeds							
Service deionization (SDI) – MΩ.cm	1 MΩ.cm minimum resistivity at exhaustion							
Reverse Osmosis (RO) – µS/cm	Recommended <30 µS/cm							
Free Chlorine		0.05 ppm max						
тос	<50 ppb max (RO feed)							
Carbon dioxide		30 ppm (max recommended)						
Silica		2 ppm (max recommended)						
Particulates	Filtration down to 5-10 mi	cron advisable to protect internal and/or po	pint of use filters					
Temperature		1-35°C (Recommend 10-15°C)						
Flowrate (maximum requirement)	130 l/hr (34 USG)							
Drain requirements	Up to 2 l/min (0.5 USG)							
Feedwater pressure	0.7 bar (10 psi) maximum; 0.07 bar (1 psi) minimum							
* Fit LA652 Pressure Regulator where feedwater pressure e	xceeds specified limits							
Dimensions	Heigh	435mm, Width 375mm, Depth 340mm						
Weight	19 kg (42 lbs) 19 kg (42 lbs) 18 kg (40 lbs)							



Process Flow PURELAB Chorus 1 Complete



Process Flow Diagram PURELAB Chorus 1 Complete With Real-Time TOC Monitoring



Application		PURELAB Chorus Complete				
Nominal output at 15°C		10 l/hr	20 l/hr			
Dispense Flowrate		>1.5 l/min	>1.5 l/min			
Inorganics @25°C		18.2 MΩ.cm	18.2 MΩ.cm			
Total organic carbon (TOC)		<5 ppb	<5 ppb			
Bacteria		<0.001 CFU/ml°	<0.001 CFU/ml°			
Bacterial Endotoxin		<0.001 EU/ml <sup>*</sup>	<0.001 EU/ml <sup>*</sup>			
рН		Effectively neutral	Effectively neutral			
Particles (filtration)		0.2 μm <sup>1</sup>	0.2µm'			
RNase		<1 pg/ml	<1 pg/ml			
DNase		<5 pg/ml	<5 pg/ml			
Purification pack capacity		Litres to 18.2 MΩ.cm = 94,100/(μ	S/cm + (2.3 x ppm CO_))			
° With LC134/145/197 POU filter/Bi	ofilter <sup>+</sup> With LC197 Biofilter <sup>+</sup> With LC1	.34/145 POU filter				
	Inorganics @25 C	>1MQ.cm, typically >10MQ.cm				
Reservoir Specification (Type 2)*	тос	50ppb, typically <25ppb				
	Bacteria	100CFU/ml, typically <20CFU/ml				
*Typical values for reservoir water						
Source		Portable Mains Wat	er Supply			
Fouling index (max)		<10				
Free Chlorine		0.5 ppm max				
Carbon dioxide		ldeally <20 ppm				
Silica		30 ppm (max recommended)				
Temperature		1-35°C (Recommend 10-15°C)				
Flowrate (maximum requirement)		130 l/hr (34 USG)				
Drain requirements		Up to 2 1/min (0.5 USG)				
Feedwater pressure		4.0 bar (60 psi) min; 6 bar (90 psi) max* With boost pump: flooded suction (min) to 2.0 bar (30 psi) max				
* Fit LA652 Pressure Regulator when	e feedwater pressure exceeds specified	limits				
Dimensions		Height 679mm, Width 376mm, Depth 353mm				
Weight (with boost pump)		17 kg (38 lbs)	18 kg (40 lbs)			
Weight		15 kg (33 lbs)	16 kg (36 lbs)			

### Pharma Compliance

#### Pharma Compliance Flow Diagram



Application	Pharma Compliance - VCLSM1
Resistivity	18.2 MΩ.cm
Dispenser flow-rate	up to 2 L/min
Total Organic Carbon (TOC)	1 – 3 ppb
Bacterial Endotoxin	<0.001 EU/ml
Bacterial Spec	<0.001 CFU/ml with a LC134 or LC145 or LC197 Point-of-use 0.2µm Micron filter fitted.
DNase	<5 pg/ml
RNase	<1 pg/ml
Particles	<0.01 µm
Recirculation Mode	During periods of non-use the unit will automatically operate in intermittent (10 minutes every 1 hours) re-circulation mode to maintain water purity with maximum efficiency.
Water Source	Pre-treated preferably RO, SDI or distilled .
Fouling Index (max)	1 for all models. A 5 – 10 micron pre-filter is recommend for all non RO feeds.
Service Deionization (SDI)	1MΩ.cm minimum at exhaustion.
Reverse Osmosis (RO)	Recommend <30 µS/cm
Free Chlorine (max)	0.05 ppm
тос	0.05 ppm max
Carbon Dioxide	30 ppm max
Silica	2 ppm max
Particulates	Filtration down to 0.2 micron advisable.
Temperature	1 – 40°C (Recommended 10 – 15°C)
Maximum Inlet Pressure	0.7 bar (10 psi) Fit a LA652 Pressure Regulator where feedwater exceeds specified limits.

Maximum Inlet Pressure	0.7 bar (10 psi) Fit a LA652 Pressure Regulator where feedwater exceeds specified limits.
Minimum Inlet Pressure	0.07 bar (1 psi)
Flowrate	130 l/hr (34 USG)
Drain Requirements	Up to 2 l/min (0.5 USG) (Gravity fall with air gap) Max during service.
timize your water purity at the point of use	

Biofilter: Endotoxin removal (<0.001 EU/ml)

Microfilter: Particulate removal (≥0.2 μm)

TREATED WATER SPECIFICATIONS



#### Process Flow PURELAB Chorus 2+ (RO/DI/UV)

Process Flow PURELAB Chorus 2+ (RO/DI/UV)

BO module for the 201 v

Degassin Module



dule for the 20I variant only

Application	PURELAB Choru	ıs 2+ (RO/DI/UV)	PURELAB Chorus 2+ (RO/EDI/UV)						
Nominal output at 15°C	10 l/hr*	20 l/hr*	10 l/hr*	20 l/hr*					
Nominal daily output	240 l/day	480 l/day	220 l/day	220 l/day					
Inorganics @25℃	10->1	5 MΩ.cm	10->1	L5 MΩ.cm					
Organics (MW>200 Dalton)	>99% 1	rejection	>99%	rejection					
Total organic carbon (TOC)	<10	) ppb	<1	.0 ppb					
Bacteria	<0.001	CFU/ml°	<0.00	1 CFU/ml°					
рН	Effective	ly neutral	Effectiv	ely neutral					
Particles (filtration)	0.2	!µm∘	0.	2µm∘					
Purification pack capacity	Litres to 15 MΩ.cm = 74,70	00/(μS/cm + (2.3 x ppm CO,))	Litres to 15 MΩ.cm = 74,7	/00/(μ5/cm + (2.3 x ppm CO <u>,</u> ))					
* Standard conditions are 4 bar inlet pressure at 1 Refer to flow tables outside these conditions. ° Wi	* Standard conditions are 4 bar inlet pressure at 15°C, fed with potable water and a clean pre-treatment cartridge. Refer to flow tables outside these conditions. ° With LC134/145 POU filter								
Source	Potable main	is water supply	Potable mains water supply						
Fouling index (max)	<	10	<10						
Conductivity	<2000	) μS/cm	<2000 µS/cm						
Free Chlorine	0.5 pp	om max	0.5 ppm max						
Heavy Metals (max)	0.05	5 ppm	0.05 ppm						
Silica	30	ppm	30 ppm						
Temperature	1-	35℃	1-35℃						
Flowrate (maximum requirement)	100 l/h	r (27 USG)	100 l/hr (27 USG)						
Drain requirements	80 l/hr	(21 USG)	80 l/hr (21 USG)						
Feedwater pressure	4.0 bar (60 psi) mir With boost pump: flooded suc	n; 6 bar (90 psi) max* tion (min) to 2.0 bar (30 psi) max	4.0 bar (60 psi) min; 6 bar (90 psi) max* With boost pump: flooded suction (min) to 2.0 bar (30 psi) max						
*Fit LA652 regulator where feedwater pressure ex	ceeds specified limits								
Dimensions	Height 679mm, Width	376mm, Depth 353mm	Height 679mm, Widt	h 376mm, Depth 353mm					
Weight (with boost pump)	17 kg (37 lbs)	18 kg (40 lbs)	17 kg (37 lbs)	18 kg (40 lbs)					
Weight	15 kg (33 lbs)	16 kg (35 lbs)	15 kg (33 lbs) 16 kg (35 lb						



#### Process Flow PURELAB Chorus 2 (RO/DI)

Process Flow PURELAB Chorus 3 (RO)



Application	PURELAB Cho	PURELAB Chorus 3 (RO)					
Nominal output at 15°C	10 l/hr	20 l/hr	10 l/hr	20 l/hr	30 l/hr		
Nominal daily output	240 l/day	480 l/day	240 l/day	480 l/day	720 l/day		
Inorganics @25°C	1 to >10	) MΩ.cm		>95% rejection			
Organics (MW>200 Dalton)	>99% r	ejection		>95% rejection			
Total organic carbon (TOC)	<30	ррb		<50 ppb			
Bacteria*	<5 Cl	-U/ml		<50 CFU/ml			
рН	Effective	ly neutral		Effectively neutral			
Particles	>99% r	ejection		>99% rejection			
Purification pack capacity	Litres to 1MΩ.cm + (2.3 x	= 103,200/(µS/cm opm CO <u>,</u> ))		N/A			
Standard conditions are 4 bar inlet pressure at 15°C, fed with potable water and a clean pre-treatment cartridge. Refer to flow tables outside these conditions. * Bacterial specification is before the water goes into the reservoir							
Source – originally from potable supply, then pretreated	Potable main	s water supply	Potable mains water supply				
Fouling index (max)	1	10					
Conductivity	<2000	μS/cm	<2000 µS/cm				
Free Chlorine	0.5 pp	m max	0.5 ppm max				
Heavy Metals (max)	0.05	ppm	0.05 ppm				
Silica	30	ppm	30 ppm				
Temperature	1-3	35°C	1-35°C				
Flowrate (maximum requirement)	100 l/hr	(27 USG)	100 l/hr (27 USG)				
Drain requirements	80 l/hr		80 l/hr (21 USG)				
Feedwater pressure	2.0 bar (30 p 0.5 bar (7.5 p	2.0 bar (30 psi) maximum; 0.5 bar (7.5 psi) minimum**					
**Fit LA652 Regulator where feedwater pressure exceed	s specified limits						
Dimensions Height 435mm, Width 375mm, Depth 340mm							

# Find your product

						6 <b>1</b> .	6 <b>1</b> )	ß	ß	ĥ					
		C	horus	1			Choru	s 2 & 3			fle	ex		Qu	est
	9	l Research	cience	te	ompliance	(\/U/IC	(VU/I							+	Quest
	Life Scien	Analytica	General <u>5</u>	1 Comple	Pharma (	2+ (RO/EI	2+ (RO/D	2 (RO/DI)	3 (RO)	flex 1	flex 2	flex 3	flex 3+	Ŋ	Non UV
Water Type															
Ultrapure Type I	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$					√#	$\checkmark$	$\checkmark$	$\checkmark$	~	~
Pure Type II				$\checkmark$		$\checkmark$	$\checkmark$	$\checkmark$		√*		$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$
RO Type III									$\checkmark$	√*				~	~
Impurities to remove															
Nucleases	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$				$\checkmark^{\diamond}$	$\checkmark$	$\checkmark$	$\checkmark$	
Endotoxins/Pyrogens	$\checkmark$	√◊	$\checkmark^{\diamond}$	$\checkmark^{\diamond}$	$\checkmark$	$\checkmark^{\diamond}$	$\checkmark^{\diamond}$			√#		$\checkmark$ $\diamond$	$\checkmark$ $\diamond$	√ ◊	√ ◊
Inorganics	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	~	$\checkmark$
Organics	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$
Bacteria	$\checkmark$ $\land$	$\checkmark$ $\land$	$\checkmark^{\Delta}$	$\checkmark^{\Delta}$	$\checkmark$ $\land$	$\checkmark^{\Delta}$	$\checkmark$ $\land$	$\checkmark$ $\land$	$\checkmark$	$\checkmark^{\Delta}$	$\checkmark^{\Delta}$	$\checkmark$	$\checkmark$ $\land$	√∆	V۵
Particulates	$\checkmark$	$\checkmark$	$\checkmark^{\Delta}$	$\checkmark$ $\land$	$\checkmark$ $\land$	$\checkmark$ $\land$	$\checkmark$ $\land$	$\checkmark$ $\land$	$\checkmark$ $\land$	$\checkmark \land$	$\checkmark$	√ ∆	$\checkmark$	$\checkmark^{\Delta}$	√∆
Features															
PureSure®	$\checkmark$	$\checkmark$	$\checkmark$		$\checkmark$										
Real time TOC monitoring	$\checkmark$	$\checkmark$		√.	$\checkmark$						$\checkmark$	$\checkmark$	$\checkmark$		
Potable tap water				$\checkmark$		$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$			$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$
Wall mounting	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	~	~
Floor mounting	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$						
Purity monitoring to POU	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$			$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	~	~
Halo Dispense compatible	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$								
Remote Dispense Compatible	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$							$\checkmark$	$\checkmark$
Auto Volume Dispense	$\checkmark^{\dagger}$	$\checkmark^{\dagger}$	$\checkmark^{\dagger}$	$\checkmark^{\dagger}$	$\checkmark$	$\checkmark^{\dagger}$	$\checkmark^{\dagger}$			$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$
Variable flow rate dispense	$\checkmark^{\dagger}$	$\checkmark^{\dagger}$	$\checkmark^{\dagger}$	$\checkmark^{\dagger}$	$\checkmark$	$\checkmark^{\dagger}$	$\checkmark^{\dagger}$			$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$		
Drop-by-drop control	$\checkmark^{\dagger}$	$\checkmark^{\dagger}$	$\checkmark^{\dagger}$	$\checkmark^{\dagger}$	$\checkmark$	$\checkmark^{\dagger}$	$\checkmark^{\dagger}$			$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$		
Locked dispense	$\checkmark^{\dagger}$	$\checkmark^{\dagger}$	$\checkmark^{\dagger}$	$\checkmark^{\dagger}$	$\checkmark$	$\checkmark^{\dagger}$	$\checkmark^{\dagger}$			$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$
USB connection	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	~	$\checkmark$	$\checkmark$	$\checkmark$
Full product validation*	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$		

<sup>†</sup> When fitted with a Halo dispenser solution

\*\* Full product validation needs to be purchased separately.  $\checkmark \bullet$  Only PC120COBPM1-TOC

# Find your product

	7	7	7	4		1,	1,	6	41 °	ſ	l				
			Chorus	1			Choru	<b>is 2 &amp;</b> 3	3		PUREL	AB fle>	<b>(</b>	Qı	lest
	Ge	l Research	cience	ė	ompliance	(VU/)(	/UV)								Quest
	Life Scien	Analytica	General S	1 Complet	Pharma C	2+ (RO/ED	2+ (RO/DI	2 (RO/DI)	3 (RO)	flex 1	flex 2	flex 3	flex 3+	۸N	Non UV
Ideal solution for Systems also have wider applicabi	lity. Spea	ιk to yoι	ur local E	LGA spec	cialist fo	r furthe	r inform	ation.							
Cell culture	$\checkmark$			$\checkmark$							$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$
Liquid Chromatography (HPLC, UHPLC)	$\checkmark$	~		$\checkmark$							~	~	$\checkmark$	V	
Microbiological Analysis	$\checkmark$	$\checkmark$		$\checkmark$							$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	
Genetic (PCR, DNA/RNA sequencing, DNA, Nucleic acid)	$\checkmark$														
Gas Chromatography		$\checkmark$		$\checkmark$							$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	
Electrochemistry		$\checkmark$	$\checkmark$	$\checkmark$		$\checkmark$	$\checkmark$				$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	
Immunochemistry	$\checkmark$			$\checkmark$							$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	
Atomic Spectroscopy (Flame AA, GFAA, ICP-AE)		$\checkmark$	$\checkmark$	$\checkmark$							$\checkmark$	$\checkmark$	$\checkmark$	V	$\checkmark$
Mass Spectrometry (ICP-MS, GC-MS, LC-MS)	$\checkmark$	$\checkmark$		$\checkmark$							$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	
General lab water requirement (glassware washing, heating baths, autoclave filling)				$\checkmark$				V	$\checkmark$	$\checkmark$		~	V	$\checkmark$	$\checkmark$
Spectrophotometry (inc. UV, IR, near UV, nearIR)				$\checkmark$		$\checkmark$	$\checkmark$				$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	
Feed to ultrapure water system						$\checkmark$	$\checkmark$	~	$\checkmark$						
Media/buffer preparation (inc pH solution)						$\checkmark$	$\checkmark$	$\checkmark$						$\checkmark$	~
General chemistry (inc Titrimetry)			$\checkmark$	~		$\checkmark$	$\checkmark$	~	$\checkmark$	~	~	~	~	$\checkmark$	~

### **Product Part Numbers**



Part Number	
PF1-M3	PURELAB flex 1 (optional Purification Pack)
PF2-M3	PURELAB flex 2 (Purification Pack + UV + TOC)
PF3-M3	PURELAB flex 3
PF3+-M3	PURELAB flex 3+

Unit is supplied as standard with: All units: quick Reference Guide & user manual, certificate of conformity, bypass pack. Flex 2: 1 x LC210. Flex 3 & 3+: 1 x LC217, 1 x LC210, 1 x LC216

PURELAB Chorus 1							
Part Number	Description						
PC1ANRXM2	PURELAB Chorus 1 Analytic Research Ultrapure System						
PC1LSCXM2	PURELAB Chorus 1 Life Science Ultrapure System						
PC1GSCXM2	PURELAB Chorus 1 General Science Ultrapure System						
Unit is supplied as standard with: Fitted/Included: appropriate consumables (including 2 x LC232 - see list below), MANU39998/MANU39997/MANU40001 Operator Manual, GUID39864							

Unit is supplied as standard with: Fitted/Included: appropriate consumables (including 2 x LC232 - see list below), MANU39998/MANU39997/MANU40001 Operator Manual, GUID39864
Light Guide, GUID40005 Quick Reference Guide, INST40012 Quick Start Guide, LA762 basic installation kit, LC233 x 2 bypass block.

PURELAB Chorus 1 Complete	Description						
Part Number							
PC110COXXM1	PURELAB Chorus 1 Complete 10 L/hr System						
PC110COBPM1	PURELAB Chorus 1 Complete 10 L/hr System with Boost Pump						
PC120COXXM1	PURELAB Chorus 1 Complete 20 L/hr System						
PC120COBPM1	PURELAB Chorus 1 Complete 20 L/hr System with Boost Pump						
PC120COBPM1-TOC	PURELAB Chorus 1 Complete 20 L/hr System with Boost Pump and TOC Monitoring						

Unit is supplied as standard with: Fitted/Included: appropriate quantity LC240, 1 x LC241, 1 x LC275, 1 x LC210, MANU40932 Operator Manual (Non-TOC models) 1 x LC322, 1 x LC241, 1 x LC275, 1 x LC323, 1 x LC323, 1 x LC210 (PC120COBPM1-TOC) GUID39864 Light Guide, LA762 Basic installation Kit.

PURELAB Pharma Compliance	Description	Dimension
Part Number	Description	Dispenser
VCLSDM1	PURELAB Pharma Compliance	1x Halo
VCLSDM1-D1	PURELAB Pharma Compliance with 1x Dispenser	1x Halo, 1x Dispenser
VCLSDM1-D2	PURELAB Pharma Compliance with 2x Dispenser Biofilter fitted	1x Halo, 2x Dispenser
VCLSDM1-D3	PURELAB Pharma Compliance with 3x Dispenser	1x Halo, 3x Dispenser
Unit is supplied as standard with: Fitted/Included: 2 x LC232, 1 x LC210, 1 x LC151, POWE40855, MANU41381 Operator Manual, VCLSDM1 Validation Manual, LA835 Installation Kit		

PURELAB Chorus 2+ RO/EDI/UV	Description	
Part Number		
PC210EUXXM1	PURELAB Chorus 2+ RO/EDI/UV 10 l/hr System	
PC210EUBPM1	PURELAB Chorus 2+ RO/EDI/UV 10 l/hr System with Boost Pump	
PC220EUXXM1	PURELAB Chorus 2+ RO/EDI/UV 201/hr System	
PC220EUBPM1	PURELAB Chorus 2+ RO/EDI/UV 20 1/hr System with Boost Pump	
Unit is supplied as standard with: Fitted/Included: appropriate quantity LC240, 1 x LC241,1 x LC243, 1 x LC277, 1 x LC285, 1 x LC233, MANU40932 Operator Manual, GUID39864 Light Guide, LA762 Basic installation Kit		

PURELAB Chorus 2+ RO/DI/UV Part Number	Description	
PC210DUXXM1	PURELAB Chorus 2+ RO/DI/UV 10 l/hr System	
PC210DUBPM1	PURELAB Chorus 2+ RO/DI/UV 10 l/hr System with Boost Pump	
PC220DUXXM1	PURELAB Chorus 2+ RO/DI/UV 20 1/hr System	
PC220DUBPM1	PURELAB Chorus 2+ RO/DI/UV 20 1/hr System with Boost Pump	
Unit is supplied as standard with: Fitted /Included: appropriate quantity LC240, 1 x LC212, 1 x LC272, 1 x LC274, 1 x LC275, MANU40932 Operator Manual,		

GUID39864 Light Guide, LA762 Basic insallation Kit.

PURELAB Chorus 2 RO/DI	Description	
Part Number		
PC210DIXXM3	PURELAB Chorus 2 RO/DI 10 1/hr System	
PC210DIBPM3	PURELAB Chorus 2 RO/DI 10 l/hr System with Boost Pump	
PC220DIXXM3	PURELAB Chorus 2 RO/DI 20 1/hr System	
PC220DIBPM3	PURELAB Chorus 2 RO/DI 20 1/hr System with Boost Pump	
Unit is supplied as standard with: Fitted/Included: appropriate quantity LC240 RO, 1 x LC241, 1 x LC234, MANU40003, Operator Manual, GUID39864 Light Guide, GUID40005		

Quick Reference Guide, INST40009 Quick Start Guide, LA762 basic installation kit, LC233 bypass block















### Product Part Numbers

PURELAB Chorus 3	Description	
Part Number	Description	F
RO310XXM3	PURELAB Chorus 3 RO 10 l/hr	
RO310BPM3	PURELAB Chorus 3 RO 10 l/hr with Boost Pump	
RO320XXM3	PURELAB Chorus 3 RO 20 1/hr	
RO320BPM3	PURELAB Chorus 3 RO 20 1/hr with Boost Pump	
RO330XXM3	PURELAB Chorus 3 RO 30 l/hr	
RO330BPM3	PURELAB Chorus 3 RO 30 l/hr with Boost Pump	

Unit is supplied as standard with: Fitted/Included: appropriate quantity LC240 RO, 1 x LC241, GUID39864 Light Guide, GUID40004 Quick Reference Guide. INST40008 Quick Start Guide. LA762 basic installation kit

PURELAB Chorus Reservoirs	Description
Part Number	Description
LA757	15 Litre Reservoir
LA758	30 Litre Reservoir
LA759	60 Litre Reservoir
LA760	100 Litre Reservoir

Each reservoir supplied with LA773 Installation Kit & LC216 Composite Vent Filter

PURELAB Dispenser	Developing	
Part Number	Description	
LA824	PURELAB Dispenser Installation Kit	
LA827	(Worktop) Bench Mounting Bracket Kit	
LC134	Point-Of-Use 0.2µm Filter (POU). Recommended Change: 3 months	
LC145	Point-Of-Use 0.2µm Filter (POU). Recommended Change: 3 months	
LC197	Point-Of-Use Bio Filter (POU). Recommended Change: 3 months	
POWE40855	Power Supply	

Unit is supplied as standard with: 1 x LA826 Installation kit & MANU41319 Operator manual. PURELAB Dispenser Installation kit & accessories

PURELAB Chorus Halo Dispensers		
Part Number	Description	
LA754	Halo Dispenser	1
LA755	Advanced Halo Dispenser	-
LA756	Flexible Dispenser	

Each supplied with GUID40007 or GUID40006 Quick Reference Guide, INST40011 or INST40013 Quick Start Guide, LA774 Installation Kit (plus for LA756: MANU40002 Operator manual)





ELGA 🕡 VEOLIA

# Dedicated to Discovery

#### info@elgalabwater.com / www.elgalabwater.com

ELGA Labwater are specialists in the engineering, service & support of water purification systems.

Unrivalled product design has achieved international recognition and awards.

Worldwide technical service teams support science & healthcare globally with specialist expertise.

Global digital performance monitoring from Hubgrade ensures laboratory work is uninterrupted.

A global supply chain supports clients from regional centres on every continent.

To find your nearest ELGA representative, go to *www.elgalabwater.com* and select your country for contact details.

**Elga Global Operations Centre.** tel: +44 (0) 203 567 7300 fax: +44 (0) 203 567 7205













OVER 80 INTERNATIONAL PATENTS

ELGA is the global laboratory water brand name of Veolia. The information contained in this document is the property of VWS (UK) Ltd, trading as ELGA LabWater, and is supplied without liability for errors or omissions. © VWS (UK) Ltd. 2022 - All rights reserved. ELGA® and MEDICA® are registered trademarks of VWS (UK) Ltd. P0134 PURELAB Range Brochure