
# PURELAB<sup>®</sup> Chorus<sup>®</sup> 3

*Operator Manual*  
MANU39996 VER 02 3/17

## WELCOME

Thank you for purchasing the Chorus 3. Please read the following Health and Safety notice before proceeding.

### HEALTH AND SAFETY

Please follow the safety information detailed in this instruction sheet. Failure to observe these instructions could result in damage to the equipment, processes, associated equipment and injury or death to persons.

#### 1. Electricity

##### WARNING!

POSITION THE POWER SUPPLY SO THAT IT CANNOT COME INTO CONTACT WITH WATER.

#### 2. Residual Preservative

##### WARNING!

DURING THE COMMISSIONING CYCLE TRACE RESIDUAL PRESERVATIVE IS FLUSHED FROM THE SYSTEM. A MATERIAL SAFETY DATA SHEET IS AVAILABLE UPON REQUEST.


#### Control of Substances Hazardous to Health (COSHH)

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
Website: [www.elgalabwater.com](http://www.elgalabwater.com)

ELGA® is the global laboratory water brand name of Veolia. PURELAB® is an ELGA Labwater trademark and technology. Owing to a policy of continual improvement, we reserve the right to amend the specifications given in this guide.  
© 2017 ELGA® Labwater / VWS (UK) Ltd. All rights reserved.


SYSTEM CONFIGURATION


Example one: 1 x PURELAB® Chorus 3 (RO) + 15L or 30L Reservoir.

Example two: 2 x PURELAB® Chorus 3 (RO) + 60L Reservoir.

Example three: 4 x PURELAB® Chorus 3 (RO) + 100L Reservoir.

Further advice on configuration connections.


**Ia** MODULAR STACKING BENCH MOUNTED, RESERVOIR ON TOP**WARNING!**


ONLY PLACE ONE CHORUS SYSTEM ONTO TOP OR BELOW THE RESERVOIR.

**1b MODULAR STACKING, RESERVOIR UNDERNEATH**


a. Screw down the reservoir plate.


b. Place on top of reservoir, line up with bolt heads. Screw in joining bolt.

**2 CONNECT THE WATER (PORT No. 1)**

a. Remove the transit plugs from the water inlet connection.


b. Using the tube supplied, firmly push one end securely into the water connector.

**3 CONNECT THE DRAIN (PORT No.3)**

a. Remove the transit plugs from the water outlet connection.


**CAUTION!**

Ensure drain pipe is not in waste water.


b. Using the tube supplied, firmly push one end securely into the water connector.

## 4 CONNECT THE PURIFIED WATER OUTLET (PORT NO.4)

a. Remove the transit plugs from the water outlet connection.

b. Using the tube supplied, firmly push one end securely into the water connector.

**N.B This product is not suitable to operate directly with a chorus dispenser.**

## 5 COMMS CONNECTIONS


### WARNING!

ONLY USE COMMS CABLE SUPPLIED BY ELGA.

## 6 CONNECT THE POWER (NO.6 PORT)

a. Connect the power cable as shown.

b. Only use the power supply supplied with the Chorus. The power supply must be earthed.

c. **PRESS ONCE** to start Commission mode. Bottom red light will be flashing until cycle has ended.

Approx 60 minutes

### WARNING!

ONLY USE THE APPLIANCE COUPLER (MAINS LEAD) AND POWER SUPPLY PROVIDED. IF THE EQUIPMENT IS USED IN A MANNER NOT SPECIFIED BY ELGA, THE PROTECTION PROVIDED BY THE EQUIPMENT MAY BE IMPAIRED. POSITION THE POWER SUPPLY SO THAT IT CANNOT COME INTO CONTACT WITH WATER.

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PURELAB®  
Chorus

## Tablet of contents

|     |                                                   |                |
|-----|---------------------------------------------------|----------------|
| 1.  | <b>Introduction</b>                               | <b>Page 1</b>  |
| 1.1 | Health & Safety                                   | Page 1         |
| 1.2 | Product Range                                     | Page 1         |
| 1.3 | Use of this Manual                                | Page 1         |
| 1.4 | Installation                                      | Page 1         |
| 1.5 | System Configuration                              | Page 1         |
| 1.6 | Environment                                       | Page 1         |
| 1.7 | Commissioning                                     | Page 1         |
| 1.8 | Customer Support                                  | Page 1         |
| 2.  | <b>Your Guide to Chorus</b>                       | <b>Page 2</b>  |
| 3.  | <b>Maintenance</b>                                | <b>Page 3</b>  |
| 3.1 | Cleaning the inlet filter assembly                | Page 3         |
| 3.2 | Replacing the Pre-treatment pack                  | Page 4         |
| 3.3 | Replacing the RO Reverse osmosis module           | Page 5         |
| 3.4 | Replacing the Composite vent filter               | Page 6         |
| 4.  | <b>Self Help Guidelines</b>                       | Page 6         |
| 5.  | Sanitisation Procedure                            | Page 7 - 8     |
| 6.  | <b>Health &amp; Safety notes</b>                  | <b>Page 9</b>  |
| 6.1 | Environment                                       | Page 9         |
| 6.2 | Electricity                                       | Page 9         |
| 6.3 | Residual Biostat                                  | Page 9         |
| 6.4 | Control Of Substances Hazardous to Health (COSHH) | Page 9         |
| 7.  | <b>Technical Specifications</b>                   | <b>Page 10</b> |
| 8.  | <b>Customized Operation</b>                       | <b>Page 11</b> |
| 8.1 | USB Ticket Printing                               | Page 11        |
| 9.  | <b>Warranty / Conditions</b>                      | <b>Page 12</b> |

## Copyright Note

The information contained in this document is the property of VWS (UK) Ltd, trading as ELGA® LabWater, and is supplied without liability for errors or omissions.

No part of this document may be reproduced or used except as authorized by contract or other written permission from VWS (UK) Ltd. The copyright and all restrictions on reproducing and use apply to all media in which this information may be placed.

VWS (UK) Ltd. Pursue a policy of continual product improvement and reserve the right to alter without notice the specification, design, price or conditions of supply of any product or service.

© VWS (UK) Ltd. 2017 - All rights reserved.

ELGA® and PURELAB® Are registered trademarks of VWS (UK) Ltd.

## Copyright Note

Effersan™ is a registered trademark of Activon Inc.

## **INTRODUCTION**

### **1.1 Health & Safety**

Please ensure you read the Health & Safety notes in Section 6.

### **1.2 Product Range**

This Operator manual has been prepared for Chorus product models:

- PURELAB® Chorus Part No. RO310XXM2
- PURELAB® Chorus Part No. RO320XXM2
- PURELAB® Chorus Part No. RO330XXM2
- PURELAB® Chorus Part No. RO310BPM2
- PURELAB® Chorus Part No. RO320BPM2
- PURELAB® Chorus Part No. RO330BPM2

### **1.3 Use of this Manual**

This manual guides you through the operation and maintenance of the Chorus allowing you to obtain a guaranteed supply of purified water to meet your requirements.

### **1.4 Installation**

Chorus is supplied with a Quick Start Guide.

Quick Start Guide shows you how to install and commission the Chorus.  
See Quick Start Guide for typical configurations.

### **1.5 System Configuration**

Maximum system size 4 x Chorus 3 RO + Reservoir.  
See Quick Start Guide for typical configurations.

### **1.6 Environment**

The Chorus unit should be installed on a flat, level surface, in a clean, dry environment.  
The unit can also be wall mounted using the purpose designed wall mounting kit against a vertical wall capable of supporting the combined weight of the unit and mounting kit.

#### **WARNING!**

NOT TO BE USED IN FUME EXTRACTION ENCLOSURES.

### **1.7 Commissioning**


Chorus is supplied in a preset commissioning mode which must be completed before purified water can be dispensed. Refer to Self Help Guidelines, section 4 Process Button.

### **1.8 Customer Support**


If you need help with your Chorus, please call your local ELGA® representative. For the address of the nearest ELGA® LabWater Sales and Service office visit the country list on our website.  
[www.elgalabwater.com](http://www.elgalabwater.com) or contact ELGA® Lab-Water at: Email: [techsupport@elgalabwater.com](mailto:techsupport@elgalabwater.com)


2 YOUR GUIDE TO THE CHORUS

- Status Light
- Warning Light
- Process Button
- RO Modules


Use a slotted screw driver to release quarter turn slots on Chorus.

- Sanitization Port
- Pre-treatment Pack

- 9. COMMS 2
- 8. Power In 2
- 7. COMMS 1
- 6. Power In 1
- 5. Return
- 4. Product Outlet
- 3. Drain
- 1. Water Inlet

## 3 MAINTENANCE

### WARNING!

An approved supplier or distributor should carry out any maintenance work not included in this manual.

### Note

Note: Disposal of all end-of-life consumable items should be in accordance with local statutory regulations.

### 3.1 Cleaning the Inlet Filter Assembly

Follow the instructions below.


\*Not supplied

1. Isolate the electrical supply.      2. Turn off water supply, remove inlet pipe.

3. Unscrew the inlet filter and remove the mesh filter.      4. Reconnect the water and power supply. Replace or clean and re-assemble.

### 3.2 Replacing the Pre-treatment Pack LC241


Follow the instructions below.




1. Process button off.


2. Unclip Pre-treatment cartridge LC241.

3. Disconnect and remove Pre-treatment cartridge.

4. Remove transit plugs and insert new Pre-treatment cartridge.

**WARNING!**

MODULES WILL REQUIRE COMMISSIONING PRIOR TO USE.


### 3.3 Replacing the Reverse Osmosis Module LC240

Follow the instructions below.

Chorus is fitted with Reverse Osmosis (RO) modules. The RO Modules may need replacing if the permeate water purity or flow rate is not adequate and does not meet predicted or previous performance.

Colour tubes for illustration purposes only.

#### WARNING!

ALL NEW RO MODULES ARE FILLED WITH A PRESERVATIVE SOLUTION TO PREVENT BACTERIAL CONTAMINATION DURING STORAGE. THE RO MODULES WILL THEREFORE REQUIRE RINSING PRIOR TO USE.

### 3.4 Replacing the Composite Vent Filter LC216

Follow the instructions below.

## 4 SELF HELP GUIDELINES

### Display Lights

|          |             |                                                      |
|----------|-------------|------------------------------------------------------|
| Constant | White       | Signifying water purity, unit is operating correctly |
| Blinking | White       | Consumables are nearing the end of useful life       |
| Constant | Red         | System needs immediate attention                     |
| Blinking | Red         | System in clean cycle (Sanitization / Commissioning) |
| Blinking | Red / White | Software update                                      |


### Process Button Operation

|                                                |                                      |                                                                                                                                                                                                        |
|------------------------------------------------|--------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Press<br>(At any time)                         | Start/Stop Modules process           | Start and stops modules process                                                                                                                                                                        |
| Press and hold for 15 seconds<br>(At any time) | Enter / Exit Commissioning Procedure | To enter or exit commissioning mode / sanitization procedure. Press process button for fifteen seconds this will enter commissioning mode. Press once at anytime to pause / resume commissioning mode. |
| Press Once<br>(At any time)                    | Load firmware from USB Flash devices | Automatically loads system firmware from USB flash device. Re-start on completion.                                                                                                                     |

### Default Alarm Settings

|                                    |                                                                                                                                                                                                                                                                                                   |
|------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Set Temperature Alarm              | Set default values from firmware file. Adjustable set point using application service module. Service module is an engineers tool, only an engineer can adjust the settings. Select to either 30°C, 35°C, 40°C. Default 35°C.                                                                     |
| Set Water Purity Alarm<br>Chorus 3 | Set default values from firmware file. Adjustable set point using application service module. Conductivity 20 - 100Ms/cm, Default 50 Ms/cm. Default Alarm point 15.0MΩ.cm.                                                                                                                        |
| Over fill Protection Alarm         | Feature activates drain valve if water level is greater than maximum and stops water entering the system to prevent a flood. When alarm is active: V2 energised. If alarm for 5mins, stop module process. Enter refill cycle on overflow alarm deactivation. Alarm status removed (Reservoir 95%) |




## 5 SANITISATION PROCEDURE

**1.** Process button off

**2.** Change Pre-treatment cartridge LC241. (Optional)

**3.** Insert a new LC241 and dispose of the old pack.


- 4.** Add CT1 or Effersan tablet to the sanitisation port.


Chlorine Tablets:

CT1 Chlorine tablets or EfferSan™ multi-purpose disinfecting tablets are required to successfully complete this procedure.  
(United States EPA Approved EfferSan™ Tablets)

EfferSan™ are not supplied and must be sourced from your local ELGA® LabWater distributor. Refer to the EfferSan™ label for additional information and follow all applicable directions for use on the manufacturers label in conjunction with the following instructions.

- 5.** Press and hold process button for 15 seconds to enter procedure. (commissioning)

- 6.** System will rinse for 60 mins.

- 7.** After 60 mins system automatically enters normal operation.

## 6 HEALTH & SAFETY NOTES

### **WARNING!**

WARNINGS ARE GIVEN WHERE FAILURE TO OBSERVE THE INSTRUCTIONS COULD RESULT IN INJURY OR FATALITY.

### **CAUTION!**

Cautions are given where failure to observe the instructions could result in damage to the equipment, associated equipment and processes.

### 6.1 Environment

### **CAUTION!**

Clean dry indoor. Temp 5-40°C. Humidity max 80% non-condensing.  
Failure to follow the environmental specification could result in damage to the system.  
The Chorus is not designed for use in fume cupboards where chemicals could damage the system.

### **WARNING!**

ALL NEW RO MODULES ARE FILLED WITH A PRESERVATIVE SOLUTION TO PREVENT BACTERIAL CONTAMINATION DURING STORAGE.  
THE RO MODULES WILL THEREFORE REQUIRE RINSING PRIOR TO USE.

### 6.2 Electricity

The appliance coupler (mains lead ) or power supply connected to the rear of the unit can be removed to isolate the power supply. If access to this is restricted then it is recommended that access to the supply socket is easily available to disconnect the electrical supply.

### **WARNING!**

ONLY USE THE APPLIANCE COUPLER (MAINS LEAD ) AND POWER SUPPLY PROVIDED.  
THE USE OF THESE WILL ENSURE ADEQUATE EARTH PROTECTION IS PROVIDED.

IF THE EQUIPMENT IS USED IN A MANNER NOT SPECIFIED BY ELGA, THE PROTECTION PROVIDED BY THE EQUIPMENT MAY BE IMPAIRED.

POSITION THE POWER SUPPLY SO THAT IT CANNOT COME INTO CONTACT WITH WATER

### **WARNING!**

ALWAYS ENSURE THAT THE ELECTRICAL POWER SUPPLY IS ISOLATED BEFORE WORKING INSIDE CHORUS.

### 6.3 Residual Preservative

### **WARNING!**

DURING THE COMMISSIONING CYCLE TRACE RESIDUAL PRESERVATIVE IS FLUSHED FROM THE SYSTEM.  
A MATERIAL DATA SHEET IS AVAILABLE UPON REQUEST.

### 6.4 Control of Substance Hazardous to Health (COSHH)

Material safety data sheets are available upon request.

## 7 TECHNICAL SPECIFICATIONS

### Purified Water Specification

| Model                                               | Chorus 3<br>RO310BPM2   | Chorus 3<br>RO320BPM2   | Chorus 3<br>RO330BPM2   | Chorus 3<br>RO310XXM2   | Chorus 3<br>RO320XXM2   | Chorus 3<br>RO330XXM2   |
|-----------------------------------------------------|-------------------------|-------------------------|-------------------------|-------------------------|-------------------------|-------------------------|
| Make-up rate @15°C<br>(*4bar feed)                  | 10 l/hr                 | 20 l/hr                 | 30 l/h                  | 10 l/hr *               | 20 l/h *                | 30 l/h *                |
| Delivery flow rate -<br>maximum                     | 240L/24h                | 480L/24h                | 720L/24h                | 240L/24h                | 480L/24h                | 780L/24h                |
| Inorganic Rejection                                 | >95%<br>Rejection       | >95%<br>Rejection       | >95%<br>Rejection       | >95%<br>Rejection       | >95%<br>Rejection       | >95%<br>Rejection       |
| Organics (TOC) - Typical                            | <100ppb                 | <100ppb                 | <100ppb                 | <100ppb                 | <100ppb                 | <100ppb                 |
| Organics (MW >200 Dalton)                           | >99%<br>Rejection       | >99%<br>Rejection       | >99%<br>Rejection       | >99%<br>Rejection       | >99%<br>Rejection       | >99%<br>Rejection       |
| Bacteria - Typical                                  | <5 CFU/ml               | <5 CFU/ml               | <5 CFU/ml               | <5 CFU/ml               | <5 CFU/ml               | <5 CFU/ml               |
| Bacterial Endotoxin                                 | <0.25EU/ml              | <0.25EU/ml              | <0.25EU/ml              | <0.25EU/ml              | <0.25EU/ml              | <0.25EU/ml              |
| Particles                                           | >99%<br>Rejection       | >99%<br>Rejection       | >99%<br>Rejection       | >99%<br>Rejection       | >99%<br>Rejection       | >99%<br>Rejection       |
| Feed Water Specification                            |                         |                         |                         |                         |                         |                         |
| Source Quality                                      | Potable mains<br>supply | Potable mains<br>supply | Potable mains<br>supply | Potable mains<br>supply | Potable mains<br>supply | Potable mains<br>supply |
| Fouling Index-Maximum                               | 10                      | 10                      | 10                      | 10                      | 10                      | 10                      |
| TDS-Maximum                                         | 1000ppm                 | 1000ppm                 | 1000ppm                 | 1000ppm                 | 1000ppm                 | 1000ppm                 |
| Conductivity                                        | <1400µS/cm              | <1400µS/cm              | <1400µS/cm              | <1400µS/cm              | <1400µS/cm              | <1400µS/cm              |
| Free/total Chlorine-Maximum                         | 0.5ppm                  | 0.5ppm                  | 0.5ppm                  | 0.5ppm                  | 0.5ppm                  | 0.5ppm                  |
| Chloramine - Maximum                                | 0.5ppm                  | 0.5ppm                  | 0.5ppm                  | 0.5ppm                  | 0.5ppm                  | 0.5ppm                  |
| Heavy Metals-Maximum<br>(Fe/Mn)                     | 0.05ppm                 | 0.05ppm                 | 0.05ppm                 | 0.05ppm                 | 0.05ppm                 | 0.05ppm                 |
| Silica-Maximum                                      | 30ppm                   | 30ppm                   | 30ppm                   | 30ppm                   | 30ppm                   | 30ppm                   |
| Temperature                                         | 1-35°C                  | 1-35°C                  | 1-35°C                  | 1-35°C                  | 1-35°C                  | 1-35°C                  |
| Flowrate (max requirement<br>@ 15°C)                | 89l/h                   | 100l/h                  | 103/h                   | 89l/h                   | 100l/h                  | 103/h                   |
| Drain requirements                                  | 80l/h                   | 80l/h                   | 80l/h                   | 80l/h                   | 80l/h                   | 80l/h                   |
| Maximum feed pressure<br>(w/o internal Boost pump)  |                         |                         |                         | 6.0bar<br>(90psi)       | 6.0bar<br>(90psi)       | 6.0bar<br>(90psi)       |
| Minimum feed pressure<br>(w/o internal Boost pump)  |                         |                         |                         | 4.0bar<br>(60psi)       | 4.0bar<br>(60psi)       | 4.0bar<br>(60psi)       |
| Maximum feed pressure<br>(with internal Boost pump) | 2.0bar<br>(30psi)       | 2.0bar<br>(30psi)       | 2.0bar<br>(30psi)       |                         |                         |                         |
| Minimum feed pressure<br>(with internal Boost pump) | 0.5bar<br>(7.5psi)      | 0.5bar<br>(7.5psi)      | 0.5bar<br>(7.5psi)      |                         |                         |                         |

### Connections

|                         |                    |
|-------------------------|--------------------|
| Inlet                   | 8mm (5/16) OD tube |
| Reservoir Overflow      | 8mm (5/16) OD tube |
| Drain (RO Concentrate)  | 8mm (5/16) OD tube |
| Reservoir               | 8mm (5/16) OD tube |
| Outlet (RO Concentrate) | 8mm (5/16) OD tube |

## Electrical Requirements

| Model                           | Chorus 3<br>RO310BPM2   | Chorus 3<br>RO320BPM2   | Chorus 3<br>RO330BPM2   | Chorus 3<br>RO310XXM2   | Chorus 3<br>RO320XXM2   | Chorus 3<br>RO330XXM2   |
|---------------------------------|-------------------------|-------------------------|-------------------------|-------------------------|-------------------------|-------------------------|
| Mains Input                     | 100-240 Vac<br>50-60 Hz | 100-240 Vac<br>50-60 Hz | 100-240 Vac<br>50-60 Hz | 100-240 Vac<br>50-60 Hz | 100-240 Vac<br>50-60 Hz | 100-240 Vac<br>50-60 Hz |
| System control voltage          | 24V dc                  | 24V dc                  | 24V dc                  | 24V dc                  | 24V dc                  | 24V dc                  |
| Power consumption (peak demand) | 60VA                    | 60VA                    | 60VA                    | 46VA                    | 46VA                    | 46VA                    |
| Reservoir Level Connection      | RJ45                    | RJ45                    | RJ45                    | RJ45                    | RJ45                    | RJ45                    |
| Noise level                     | <45dBA                  | <45dBA                  | <45dBA                  | <45dBA                  | <45dBA                  | <45dBA                  |

## Dimensions & Weight

| Model            | Chorus 3<br>RO310BPM2 | Chorus 3<br>RO320BPM2 | Chorus 3<br>RO330BPM2 | Chorus 3<br>RO310XXM2 | Chorus 3<br>RO320XXM2 | Chorus 3<br>RO330XXM2 |
|------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|
| Height           | 435mm                 | 435mm                 | 435mm                 | 435mm                 | 435mm                 | 435mm                 |
| Width            | 375mm                 | 375mm                 | 375mm                 | 375mm                 | 375mm                 | 375mm                 |
| Depth            | 340mm                 | 340mm                 | 340mm                 | 340mm                 | 340mm                 | 340mm                 |
| Supply Weight    | 17kg(37lb)            | 18kg(40lb)            | 19kg(42lb)            | 15kg(33lb)            | 16kg(35lb)            | 17kg(37lb)            |
| Operating Weight | 18kg(40lb)            | 19kg(42lb)            | 20kg(44lb)            | 16kg(35lb)            | 17kg(37lb)            | 18kg(40lb)            |
| Installation     | Bench/Wall/<br>Floor  | Bench/Wall/<br>Floor  | Bench/Wall/<br>Floor  | Bench/Wall/<br>Floor  | Bench/Wall/<br>Floor  | Bench/Wall/<br>Floor  |

Dependent on feed water as part of our policy of continual improvement we reserve the right to alter the specifications given in this document.

## Consumables

| Part No.  | Description           | Typical Service Life * | Max. Shelf Life |
|-----------|-----------------------|------------------------|-----------------|
| LC241     | Pre-treatment Pack    | 6 Months               | 2 years         |
| LC216     | Composite Vent Filter | 6 Months               | 2 years         |
| Effersan™ | Tablets               | >8 Months              | >8 Months       |
| CT1       | Tablets               | 1 per 6 Months         | 2 years         |

\*Service life is an estimate only, and will depend on the application and feed water quality.  
Please take care to ensure your order is the correct consumable items.

## 8 CUSTOMIZED OPERATION

It's possible to customize your Chorus by downloading additional features and changing your current settings. To find out more, please visit [www.elgalabwater.com/customize](http://www.elgalabwater.com/customize)

**8.1 USB Data Logging** To download technical data, please follow instructions.

**CAUTION!**

USE A BLANK USB.

1. Insert blank USB and press the process button.
2. To stop recording data remove USB.

ELGA® LabWater is a trading name of VWS (UK) Ltd. General Limited Warranty

VWS (UK) Ltd. warrants the products manufactured by it against defects in materials and workmanship when used in accordance with applicable instructions for a period of one year from the date of shipment for the products. VWS (UK) Ltd. MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

The warranty provided herein and the data, specifications and descriptions of the VWS (UK) Ltd. products appearing in VWS (UK) Ltd. published catalogues and product literature may not be altered except by express written agreement signed by an officer of VWS (UK) Ltd. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and, if given, should not be relied upon.

In the event of a breach of the foregoing warranty, VWS (UK) Ltd. sole obligation shall be to repair or replace, at its option, any product or part thereof that proves to be defective in materials or workmanship within the warranty period, provided the customer notifies VWS (UK) Ltd. promptly of any such defect. The exclusive remedy provided herein shall not be deemed to have failed of its essential purpose so long as VWS (UK) Ltd. is willing and able to repair or replace any nonconforming VWS (UK) Ltd. product or part. VWS (UK) shall not be liable for consequential, incidental, special or any other indirect damages resulting from economic loss or property damage sustained by any customer from the use of its products.

#### VWS (UK) Ltd. Warranty

VWS (UK) Ltd. warrants the water systems manufactured by it, BUT EXCLUDING MEMBRANES AND PURIFICATION PACKS, against defects in materials and workmanship when used in accordance with the applicable instructions and within the operating conditions specified for the systems for a period of one year from the earlier of:

the date of installation, or the 120th day following the date of shipment.

VWS (UK) LTD. MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. The warranty provided herein and the data, specifications and descriptions of the VWS (UK) Ltd. systems appearing in VWS (UK) Ltd. published catalogues and product literature may not be altered except by express written agreement signed by an officer VWS (UK) Ltd. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorised and, if given, should not be relied upon.

In the event of a breach of the foregoing warranty, VWS (UK) Ltd. sole obligation shall be to repair or replace, at its option, any product or part thereof that proves to be defective in materials or workmanship within the warranty period, provided the customer notifies VWS (UK) Ltd. promptly of any such defect. The cost of labour for the first ninety (90) days of the above warranty period is included in the warranty; thereafter, labour cost shall be at the customer's expense. The exclusive remedy provided herein shall not be deemed to have failed of its essential purpose so long as VWS (UK) Ltd. is willing and able to repair or replace any nonconforming VWS (UK) Ltd. system or component part. VWS (UK) Ltd. shall not be liable for consequential, incidental, special or any other indirect damages resulting from economic loss or property damage sustained by any customer from the use of its process systems.

Products or components manufactured by companies other than VWS (UK) Ltd. or its affiliates

("Non - VWS (UK) Ltd. products") are covered by the warranty, if any, extended by the Product manufacturer. VWS (UK) Ltd. hereby assigns to the purchaser any such warranty; however VWS (UK) LTD. EXPRESSLY DISCLAIMS ANY WARRANTY WHETHER EXPRESSED OR IMPLIED, THAT THE NON - VWS (UK) Ltd. PRODUCTS ARE MERCHANTABILITY OR FIT FOR A PARTICULAR PURPOSE.

#### NOTICE

VWS (UK) Ltd. is constantly striving to improve its products and services. Consequently, the information in this document is subject to change without notice and should not be construed as a commitment by VWS (UK) Ltd.

Also VWS (UK) Ltd. assumes no responsibility for any errors that may appear in this document. This manual is believed to be complete and accurate at the time of publication. In no event shall VWS (UK) Ltd. be liable for incidental or consequential damages in connection with or arising from the use of this manual.

VWS (UK) Ltd. warrants its products against defects in materials and workmanship as described in the Warranty statement on the preceding pages.

# The Labwater Specialists

ELGA is the global laboratory water brand name of Veolia. The information contained in this document is the property of VWS (UK) Ltd, trading as ELGA LabWater, and is supplied without liability for errors or omissions. © VWS (UK) Ltd.  
2017 - All right reserved. ELGA®, PURELAB®, ADEPT® Labpure® and PureSure® are registered trademarks of VWS (UK) Ltd.