

Think. Inside. The. Box.

PURELAB® Chorus 1 ~ General Science*



LC232/LC244
LC245/LC246

DI Pack

Lifetime expectancy: 12 months
What is purified from the water?

Ions

Impact of not changing:

Weakly bound ions are released back into the water reducing resistivity



Point of Use Filter**

LC145

Lifetime expectancy: 3 months
What is purified from the water?

Particulates and bacteria

Impact of not changing:



Potential recontamination of purified water and reduced flow rates



*A combination of treatments and technologies produce the final water specification

**For DNA, RNA & endotoxin specification at point of use please use a LC197 Biofilter

For more information please contact your local ELGA water specialist:

 +44 020 3567 7300
 info@elgalabwater.com

 e-shop!

