


ELGA  VEOLIA	EPA SPECIFICATION SHEET	EPA/E/1116-004
		REVISION 01

The following document has now been written and approved. Information relating to the document supply format is stated below. This document is intended to be circulated within ELGA ONLY.

SECTION 1		FRONT COVER	
Part Number	MANU40608	Version: 01	
Title	MR200DSM1-115-US MEDICA R200 (US only)		
Manual type	Quick Start Manual	☐	Operator Manual ☒



QRD1116	ISSUE 02	Page 1 of 3
Printed date: 31/08/16	UNCONTROLLED CHECK ISSUE BEFORE USE	

SECTION 2 EXTERNAL CARTON LABELING

SECTION 2-1 Generic position of carton labels & Carton label details

Carton Product Label



Carton UK EPA
Establishment
Number Label

**EPA Est. 091083-GBR-001
MEDICA**

ELGA 

Model MR200DSM1-115-US



Date of Manufacture 06/11/2014

Serial No. MDR00006725



Maximum Pressure

Inlet 4/60 (bar / psi)

Working 11/160 (bar / psi)

Electrical 115Vac-60Hz-2000VA

Supply



REFER TO OPERATING INSTRUCTIONS
SU REFERER AU MANUEL D'OPERATION
IN DER BEDIENUNGSANLEITUNG NACHSEHEN

ELGA® is the global laboratory water brand of
Veolia Water Solutions & Technologies.
Tel: +44 1494 887500, Web: www.elgalabwater.com
MADE IN UK

SECTION 3

PRODUCT LABEL DETAIL

SECTION 3-1


Generic position of units labels & Unit label details.

Unit Product Label



Unit UK EPA Establishment Number Label

EPA Est. No. 091083- GBR-001



AUTHORISATION - RESEARCH AND DEVELOPMENT DEPARTMENT

This document is released on behalf of ELGA and remains the property of VWS (UK) LTD.

Compiled By:

Date:

Approved:

Release Date:

QRD1116

ISSUE 02

Page 3 of 3

Printed date: 31/08/16

UNCONTROLLED CHECK ISSUE BEFORE USE