


	EPA SPECIFICATION SHEET	EPA/E/1116-023
		REVISION 02

The following document has now been written and approved. Information relating to the document supply format is stated below. This document is intended to be circulated within ELGA ONLY.

SECTION 1		FRONT COVER	
Part Number	MANU40591	Version:02	
Title	MP060RBM1-115-US – MEDICA Pro 60 (US only)		
Manual type	Quick Start Manual	☒	Operator Manual ☐



QRD1116	ISSUE 02	Page 1 of 3
Printed date: 28/11/16	UNCONTROLLED CHECK ISSUE BEFORE USE	

SECTION 2

EXTERNAL CARTON LABELING

SECTION 2-1

Generic position of carton labels & Carton label details

Carton Product Label



Carton UK EPA
Establishment
Number Label

**EPA Est. 091083-GBR-001
MEDICA**

ELGA



Model **MP060RBM1-115**



Date of Manufacture **03/09/2014**

Serial No. **MP00006667**



Maximum Pressure

Inlet **3/50** (bar / psi)

Working **11/160** (bar / psi)

Electrical Supply **115Vac-60Hz-650VA**



**REFER TO OPERATING INSTRUCTIONS
SU REFERER AU MANUEL D'OPERATION
IN DER BEDIENTUNGSANLEITUNG NACHSEHEN**

**ELGA® is the global laboratory water brand of
Veolia Water Solutions & Technologies.
Tel: +44 1494 887500, Web: www.elgalabwater.com
MADE IN UK**

SECTION 3
PRODUCT LABEL DETAIL
SECTION 3-1
Generic position of units labels & Unit label details.

 Unit
Product
Label

 Unit UK EPA
Establishment
Number Label

EPA Est. No. 091083- GBR-001

AUTHORISATION - RESEARCH AND DEVELOPMENT DEPARTMENT

This document is released on behalf of ELGA and remains the property of VWS (UK) LTD.

Compiled By:

Date:

Approved:

Release Date: