


	EPA SPECIFICATION SHEET	EPA/E/1116-054
		REVISION 02

The following document has now been written and approved. Information relating to the document supply format is stated below. This document is intended to be circulated within ELGA ONLY.

SECTION 1		FRONT COVER	
Part Number	MANU40683	Version:02	
Title	MP120RBM1-115-US – MEDICA Pro 120 (US only)		
Manual type	Quick Start Manual	☒	Operator Manual ☐



QRD1115	ISSUE 02	Page 1 of 3
Printed date: 25/06/19	UNCONTROLLED CHECK ISSUE BEFORE USE	

SECTION 2
EXTERNAL CARTON LABELING
SECTION 2-1
Generic position of carton labels & Carton label details

Carton Product Label


 Carton UK EPA
Establishment
Number Label

EPA Est. 091083-GBR-001
MEDICA

ELGA


Model MP120RXM1-115-US



Date of Manufacture 06/02/2015

Serial No. MP00006728



Maximum Pressure

Inlet 3/50 (bar / psi)

Working 11/160 (bar / psi)

 Electrical 115Vac-60Hz-650VA
Supply


REFER TO OPERATING INSTRUCTIONS
SE REFERER AU MANUEL D'OPERATION
IN DER BEDIENUNGSANLEITUNG NACHSEHEN

ELGA® is the global laboratory water brand of
Veolia Water Solutions & Technologies.
Tel: +44 1494 887500, Web: www.elgalabwater.com
MADE IN UK

SECTION 3

PRODUCT LABEL DETAIL

SECTION 3-1


Generic position of units labels & Unit label details.

Unit
Product
Label



Unit UK EPA
Establishment
Number Label

EPA Est. No. 091083- GBR-001



AUTHORISATION - RESEARCH AND DEVELOPMENT DEPARTMENT

This document is released on behalf of ELGA and remains the property of VWS (UK) LTD.

Compiled By:

Jeff Paine

Date:

25th June 2019

Approved:

[Signature]

Release Date:

25th JUNE 2019