

ELGA  VEOLIA	EPA SPECIFICATION SHEET	EPA/E/1116-077
		REVISION 01

The following document has now been written and approved. Information relating to the document supply format is stated below. This document is intended to be circulated within ELGA ONLY.

SECTION 1 FRONT COVER			
Part Number	MANU40932	Version: 01	
Title	PC220DUBPM1 – PURELAB Chorus 2+ DI 20 BP		
Manual type	Quick Start Manual ☐	Operator Manual ☒	



QRD1116	ISSUE 02	Page 1 of 3
Printed date: 10/05/17	UNCONTROLLED CHECK ISSUE BEFORE USE	


SECTION 2

EXTERNAL CARTON LABELING

SECTION 2-1

Generic position of carton labels & Carton label details

Carton UK EPA
Establishment
Number Label



Carton Product Label

**EPA Est. 091083-GBR-001
PURELAB**

ELGA



Model PC220DUBPM1



Date of Manufacture 08/05/2017

Serial No. CRS00006566



Maximum Pressure

Inlet 2.0/30 (bar / psi)

Working 6.0/90 (bar / psi)

Electrical Supply 24Vdc-155VA



REFER TO OPERATING INSTRUCTIONS
SE REFERER AU MANUEL D'OPERATION
IN DER BEDIENTUNGSANLEITUNG NACHSEHEN

ELGA® is the global laboratory water brand of
Veolia Water Solutions & Technologies.
Tel: +44 203 567 7300, Web: www.elgalabwater.com
MADE IN UK

SECTION 3

PRODUCT LABEL DETAIL

SECTION 3-1

Generic position of units labels & Unit label details.

Unit Product Label



Unit UK EPA Establishment Number Label

EPA Est. No. 091083- GBR-001

ELGA		CE	
Model PC220DUBPM1			
			
Date of Manufacture		08/05/2017	
Serial No. CRS00006566			
Maximum Pressure			
Inlet	2.0/30	(bar / psi)	
Working	6.0/90	(bar / psi)	
Electrical Supply	24Vdc-155VA		
 REFER TO OPERATING INSTRUCTIONS SE REFERER AU MANUEL D'OPERATION IN DER BEDIENTUNGSANLEITUNG NACHSEHEN			
ELGA® is the global laboratory water brand of Veolia Water Solutions & Technologies. Tel: +44 203 567 7300, Web: www.elgalabwater.com MADE IN UK			

AUTHORISATION - RESEARCH AND DEVELOPMENT DEPARTMENT

This document is released on behalf of ELGA and remains the property of VWS (UK) LTD.

Compiled By:

Date:

Approved:

Release Date:

QRD1116

ISSUE 02

Page 3 of 3

Printed date: 10/05/17

UNCONTROLLED CHECK ISSUE BEFORE USE