




| ISS | MN/PRN/DN | ACTION | APPROVED | DATE |
|-----|-----------|--------|----------|--------|
| 1 | ERN/E/422 | MBL | NSN | JAN 21 |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

| ELECTRICAL REQUIREMENTS | |
|-------------------------|---------------------|
| MAINS INPUT | 100-240 VAC 50-60Hz |
| SYSTEM VOLTAGE | 24DC |
| POWER CONSUMPTION | 30VA |
| NOISE LEVEL | <40dBA |
| UNIT WEIGHT 5Kgs | |



THIS DRAWING HAS BEEN CREATED USING CREO PARAMETRIC SOFTWARE. THE MASTER FILE IS STORED WITHIN THE WINDCHILL PDMLINK DATABASE.

IF IN DOUBT - ASK
UNCONTROLLED
Please check issue before use



MATERIAL: -
FINISH: -

GENERAL NOTE: REMOVE ALL BURRS & SHARP EDGES

TOLERANCES UNLESS STATED:
ONE PLACE DEC. ± 0.1
TWO PLACE DEC. ± 0.01
ANGULAR ± 0.5°
ALL HOLE TOLS. TO BS.4500 H11

| DIMS. IN | SCALE | DRAWN | MBL | 23/10/2020 | DATE |
|----------|-------|----------|-----|------------|------|
| mm | 1:20 | CHECKED | NSN | 15/01/2021 | |
| | | APPROVED | NSN | 15/01/2021 | |

TITLE:
ENGINEERING SALES DRAWING
LA824 STANDALONE DISPENSER

PART NUMBER: LA824

ELGA VEOLIA

ELGA LABWATER
Lane End Industrial Park
High Wycombe
Bucks, HP14 3BY, UK
ELGA LabWater is a brand name of VWS (UK) Ltd

Tel: +44 (0) 203 567 7300
Fax: +44 (0) 203 567 7305
service@elgalabwater.com

| CONNECTION DETAILS | | |
|--------------------|--------------|--------------------------------|
| 1 | INLET | 8mm (5/16") OD TUBE CONNECTION |
| 2 | POWER SOCKET | DIN SOCKET |
| 3 | POWER SOCKET | DIN SOCKET |
| 4 | COMMS | RJ45 SOCKET |
| 5 | RETURN | 8mm (5/16") OD TUBE CONNECTION |
| 6 | COMMS | RJ45 SOCKET |



THIS DRAWING IS THE PROPERTY OF ELGA LabWater. THE DETAILS SHOWN ARE STRICTLY CONFIDENTIAL AND MUST NOT BE DISCLOSED TO ANY THIRD PARTY WITHOUT THE CONSENT OF VWS (UK) LTD.

| | |
|-----------------------------|-------------------|
| DRG. NO: SD_LA824 | ISSUE 1 |
| SHEET No. 1 OF 1 | |