



ISS.	MN/PRN/DN	ACTION	APPR'D	DATE
1	DN/E/144	JP		
-	-	-	-	-
-	-	-	-	-

DOCKING VESSEL DETAILS - UNIT COMPATABILITY


UNITS SUITABLE FOR DOCKING WITH THE FOLLOWING UNITS	LA620 (35 LITRE) DOCKING VESSEL	LA621 (25 LITRE) DOCKING VESSEL
PURELAB PRIMA 7/15/20	✓	
PURELAB OPTION-S 7/15	✓	
PURELAB OPTION-R 7/15		✓
PURELAB PULSE 1 & 2		✓
MEDICA 7/15		✓
MEDICA-R 7/15		✓
MEDICA-D 7/15		✓
PURELAB ULTRA RANGE	✓	



PURELAB Ultra Range
LA620 - 35Ltr Docking Vessel.



PURELAB Prima 7, 15 or 20.
LA620 - 35Ltr Docking Vessel.



PURELAB OPTION-R 7/15.
LA621 - 25Ltr Docking Vessel.

IF IN DOUBT-ASK

3rd ANGLE PROJECTION DO NOT SCALE

MATERIAL: N/A

FINISH: N/A

GENERAL NOTE: N/A

TOLERANCES UNLESS STATED:
 ONE PLACE DEC. ± N/A ANGULAR ± 0.5°
 TWO PLACE DEC. ± ALL HOLE TOLS. TO BS.4500 H11

DIMS. IN mm	SCALE	DRAWN	JP	AUG 11	DATE
	NTS	CHECKED	SB	AUG 11	
SIZE: A3		APPROVED	LU	AUG 11	

TITLE
**SALES DIAGRAM
 TYPICAL INSTALLAT'S
 LA620/LA621
 DOCKING VESSELS**

PART NUMBER: N/A

ELGA
 www.elgalabwater.com
 info@elgalabwater.com

ELGA is the laboratory water brand of Veolia Water Solutions & Technologies.

DRG. NO. SD_DV_03	ISSUE 1
SHEET No.3 OF 3	

THIS DRAWING IS THE PROPERTY OF ELGA Labwater. THE DETAILS SHOWN ARE STRICTLY CONFIDENTIAL AND MUST NOT BE DISCLOSED TO ANY THIRD PARTY WITHOUT THE CONSENT OF VWS (UK) LTD.