

ISS.	MN/PRN/DN	ACTION	APPR'D	DATE
1	DN/E/144	JP		
-	-	-	-	-
-	-	-	-	-

DOCKING VESSEL DETAILS - ALL UNIT COMPATABILITY

UNITS SUITABLE FOR DOCKING WITH THE FOLLOWING UNITS	LA620 (35 LITRE) DOCKING VESSEL	LA621 (25 LITRE) DOCKING VESSEL
PURELAB PRIMA 7/15/20	✓	
PURELAB OPTION-S 7/15	✓	
PURELAB OPTION-R 7/15		✓
PURELAB PULSE 1 & 2		✓
MEDICA 7/15		✓
MEDICA-R 7/15		✓
MEDICA-D 7/15		✓
PURELAB ULTRA RANGE	✓	

CONNECTION & DIMENSION DETAILS.

REF	FUNCTION	DETAILS	LA621 25 LITRE
A	INLET	8mm OD TUBE CONNECTOR	✓
B	OUTLET	8mm OD TUBE CONNECTOR	✓
C	OVERFLOW	8mm OD TUBE CONNECTOR	✓
D	DISPENSE TAP	—	-
E	LEVEL CONTROL	JACK PLUG 3.5mm	✓
		JACK SOCKET 3.5mm	✓
		JACK SOCKET 1/4"	✓
F	RESERVOIR ACCESS	40mm DIAMETER	✓

UNIT WEIGHTS - ALL VERSIONS

UNIT	WEIGHT (EMPTY)	WEIGHT (FULL)
LA620 (35 LITRE)	7.0KG (15.4LB)	42.0KG (94.4LB)
LA621 (25 LITRE)	6.0KG (13.2LB)	31.0KG (68.2LB)

LEVEL CONNECTION DETAILS - ALL UNITS / ALL DOCKING VESSELS

DETAILS	JACK PLUG 3.5 mm	JACK SOCKET 3.5mm	JACK SOCKET 1/4"
PURELAB PRIMA 7/15/20	✓		
PURELAB OPTION-S 7/15	✓		
PURELAB OPTION-R 7/15	✓		
PURELAB PULSE 1 & 2	✓		
MEDICA 7/15	✓		
MEDICA-R 7/15	✓		
MEDICA-D 7/15	✓		
PURELAB ULTRA RANGE	✓	✓	

OPTIONAL WALL MOUNTING KIT - ELGA PART	No. LA591	
OVERALL HEIGHT	584mm	23.00"
DIMENSION 'A'	515mm	20.30"
DIMENSION 'B'	362mm	14.25"
DIMENSION 'C'	545mm	21.46"

IF IN DOUBT-ASK

3rd ANGLE PROJECTION DO NOT SCALE

MATERIAL: N/A




FINISH: N/A

GENERAL NOTE: N/A

TOLERANCES UNLESS STATED:
 ONE PLACE DEC. ± N/A ANGULAR ± 0.5°
 TWO PLACE DEC. ± ALL HOLE TOLS. TO BS.4500 H11

DIMS. IN mm	SCALE	DRAWN	JP	AUG 11	DATE
SIZE: A3	NTS	CHECKED	SB	AUG 11	
		APPROVED	LU	AUG 11	

LA621 - 25 LITRE DOCKING VESSEL - (DV25)



TITLE

**SALES DIAGRAM
OUTLINE & CONN'N
LA620/LA621
DOCKING VESSELS**

PART NUMBER: N/A

ELGA

www.elgalabwater.com
info@elgalabwater.com

ELGA is the laboratory water brand of Veolia Water Solutions & Technologies.

DRG. NO. **SD_DV_02** ISSUE **1**

SHEET No2 OF 3

THIS DRAWING IS THE PROPERTY OF ELGA Labwater. THE DETAILS SHOWN ARE STRICTLY CONFIDENTIAL AND MUST NOT BE DISCLOSED TO ANY THIRD PARTY WITHOUT THE CONSENT OF VWS (UK) LTD.