




ISS.	MIN/PRN/IN	ACTION	APPRD.	DATE
1	ONE/144		JP	
.
.





LA610
WALL MOUNTING KIT
(UK VERSION)




LA622
WALL MOUNTING KIT
(US VERSION)




PURELAB CLASSIC WITH LA620 DOCKING VESSEL



LA644 - INTEGRAL DISPENSE GUN
MAXIMUM REACH OF DISPENSE GUN 1 METRE



LA645 - REMOTE CONTROL STATION



LA618 - RS232 PRINTER KIT

IF IN DOUBT-ASK

3rd ANGLE PROJECTION DO NOT SCALE

MATERIAL: N/A

FINISH: N/A

GENERAL NOTE: N/A

TOLERANCES UNLESS STATED:
ONE PLACE DEC. ± N/A ANGULAR ± 0.5°
TWO PLACE DEC. ± ALL HOLE TOLS. TO BS.4500 H11

DIMS. IN	SCALE	DRAWN	JP	AUG 11
mm	NTS	CHECKED	SB	AUG 11
SIZE: A1		APPROVED	LU	AUG 11

TITLE: SALES DIAGRAM
ACCESSORIES FOR
PURELAB CLASSIC
RANGE

PART NUMBER: N/A

ELGA

www.elgalabwater.com
info@elgalabwater.com

ELGA is the laboratory water brand of Veolia Water Solutions & Technologies.

DRG. NO. SD_CL_003 ISSUE 1
SHEET No.3 OF 3

THIS DRAWING IS THE PROPERTY OF ELGA Labwater. THE DETAILS SHOWN ARE STRICTLY CONFIDENTIAL AND MUST NOT BE DISCLOSED TO ANY THIRD PARTY WITHOUT THE CONSENT OF WWS (UK) LTD.