




ISS.	MIN/PRN/DN	ACTION	APPRD.	DATE
1	DN/E/144	JP		




LA610
WALL MOUNTING KIT
(UK VERSION)




LA622
WALL MOUNTING KIT
(US VERSION)




PURELAB ULTRA WITH LA620 DOCKING VESSEL




LA618 - RS232 PRINTER KIT
100-240Vdc 50/60 Hz




LA646 - OPEN LOOP RINGMAIN KIT
MAXIMUM LENGTH OF RINGMAIN 20 METRE




LA647 - CLOSED LOOP RINGMAIN KIT
MAXIMUM LENGTH OF RINGMAIN 20 METRE



LA645 - REMOTE CONTROL STATION



LA644 - INTEGRAL DISPENSE GUN
MAXIMUM REACH OF DISPENSE GUN 1 METRE



LA643 - REMOTE DISPENSE STATION
MAXIMUM REACH OF DISPENSE GUN 2 METRE

IF IN DOUBT-ASK				
3rd ANGLE PROJECTION DO NOT SCALE				
MATERIAL: N/A				
FINISH: N/A				
GENERAL NOTE: N/A				
TOLERANCES UNLESS STATED: ONE PLACE DEC. ± N/A ANGULAR ± 0.5° TWO PLACE DEC. ± ALL HOLE TOLS. TO BS.4500 H11				
DIMS. IN	SCALE	DRAWN	JP	AUG 11
mm	NTS	CHECKED	SB	AUG 11
SIZE: A1		APPROVED	LU	AUG 11

TITLE
**SALES DIAGRAM
ACCESSORIES FOR
PURELAB ULTRA
RANGE**

PART NUMBER: N/A
ELGA
www.elgalabwater.com
info@elgalabwater.com
ELGA is the laboratory water brand of Veolia Water Solutions & Technologies.